CerTest Rapid Test

CerTest
BIOTEC

CerTest

***C. difficile* GDH+Toxin A+B**

ONE STEP *C. difficile* antigen GDH, Toxin A
and Toxin B COMBO CARD TEST

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These instructions for use apply to the following references:

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GX872001V; GX872001VC; GX872001VCC

Initial part of the Instructions for use is common for all the references/

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Specific parts for each of the references are described in the corresponding annexes (see tables below)/

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




ENGLISH

2. Intended use.....	8
3. Summary and Explanation.....	8
4. Test principle.....	9
5. Storage and Stability.....	10
6. Precautions for users.....	10
7. Specimen collection and preparation.....	11
8. Limitations of the test.....	12
9. Expected values.....	12
10. Performance characteristics.....	13
A. Analytical sensitivity (detection limit).....	13
B. Clinical sensitivity and specificity.....	13
C. Hook Effect.....	16
D. Cross Reactivity.....	16
E. Interferences.....	18
F. Repeatability and Reproducibility.....	21
ANNEX 1.....	22
I. Materials.....	22
II. Test procedures (samples).....	22
III. Interpretation of results.....	23
IV. Quality control.....	26
ANNEX 2.....	27
I. Materials.....	27
II. Test procedures (samples).....	28
III. Interpretation of results.....	29
IV. Quality control.....	31
A. External positive controls.....	32
1. CerTest C. <i>difficile</i> antigen GDH Positive Control.....	32
1.1. Positive control procedure.....	32

2. CerTest C. difficile Toxin A Positive Control	33
2.1. Positive control procedure	33
3. CerTest C. difficile Toxin B Positive Control	34
3.1. Positive control procedure	34
A.1 Interpretation positive control procedure	35
ANNEX 3	37
I. Materials	37
II. Test procedures (samples)	38
III. Interpretation of results	39
IV. Quality control	41
A. External positive controls	42
1. CerTest C. difficile antigen GDH Positive Control	42
1.1. Positive control procedure	42
2. CerTest C. difficile Toxin A Positive Control	43
2.1. Positive control procedure	43
3. CerTest C. difficile Toxin B Positive Control	44
3.1. Positive control procedure	44
A.1 Interpretation positive control procedure	45
B. External negative control. CerTest Negative Control	46
B.1 Negative control procedure	47
B.1.1 Interpretation negative control procedure	48

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ENGLISH

2. Intended use

CerTest C. difficile GDH+Toxin A+B one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH), Toxin A and Toxin B in stool samples.

CerTest C. difficile GDH+Toxin A+B combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Clostridium difficile* infection.

3. Summary and Explanation

Clostridium difficile (*C. difficile*), a Gram-positive spore bearing anaerobic bacterium is the major aetiological agent of diarrhoea and colitis associated with antibiotics (1).

C. difficile is the most common cause of health care-associated diarrhoea in developed countries and is a major source of nosocomial morbidity and mortality worldwide. Although traditionally considered an infection affecting patients in health care facilities, in the last 10-20 years *C. difficile* infection is increasingly identified in the community among people with no recent hospitalizations and without prior antibiotic exposure (2–4).

Disease due to *C. difficile* develops when the bacteria is allowed to proliferate in the colon, most commonly after antibiotic use has resulted in sufficient disturbance of the intestinal microflora. Overall, the use of antibiotics is associated with a three-fold increased risk of community-associated *C. difficile* infection (CA-CDI). However, the risk of CA-CDI is largely dependent on the class of antibiotic use (1,3).

C. difficile can release two high-molecular-weight toxins, Toxin A and Toxin B, which are responsible for the clinical manifestations, which range from mild, self-limited watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon, and death. *Clostridium difficile* Glutamate Dehydrogenase (GDH) is an enzyme produced in large quantities by all toxigenic and non-toxigenic strains, making it an excellent marker for the microorganism (2,5).

The use of a direct faecal GDH screen, together with a faecal Toxin AB test, could improve the diagnosis of *Clostridium difficile* infection. Patients with a positive faecal GDH but negative toxin can be monitored closely while awaiting the result of a toxigenic culture (TC), that is used as the gold standard technique to determine *Clostridium difficile* infection (5).

Preventing *C. difficile* infection offers a potentially significant improvement in patient's outcomes, as well as a reduction in hospital costs and resource expenditures. Careful use of antibiotics remains a prominent public health strategy to reduce *C. difficile* infection cases (1,3).

4. Test principle

CerTest C. difficile GDH+Toxin A+B is based on the principle of a qualitative immunochromatographic assay for the qualitative detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH), Toxin A and Toxin B in human stool samples from patients suspected of *Clostridium difficile* infection.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against GDH and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-GDH) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against Toxin A and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-Toxin A) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip C consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against Toxin B and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-Toxin B) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

If the sample is GDH positive, the antigen of the diluted sample reacts with the red-coloured conjugates complex (anti-GDH monoclonal antibodies-red polystyrene microspheres) in the strip A, if the sample is Toxin A positive, the antigens of the diluted sample react with the red-coloured conjugates complex (anti-Toxin A monoclonal antibodies-red polystyrene microspheres) in the strip B and if the sample is Toxin B positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-Toxin B monoclonal antibodies-red polystyrene microspheres) in strip C, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-GDH antibodies present on the membrane of strip A (test line), the anti-Toxin A antibodies present on the membrane of strip B (test line) and the anti-Toxin B antibodies present on the membrane of strip C (test line) capture the coloured conjugate and the red lines will be visible in the three strips. These bands are used to interpret the result.

If the sample is negative, there is no GDH, Toxin A and Toxin B presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-GDH, the anti-Toxin A and the anti-Toxin B antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in the three strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on three membranes will capture control green-conjugate complex and three control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

5. Storage and Stability

The kits should be stored at 2-30°C, in these conditions can be used until the expiration date which is stated on the kit label. All the kit components are for only one use and should remain in its primary packaging until use. Do not freeze.

6. Precautions for users

- **For professional *in vitro* diagnostic use only.**
- Read the instructions for use carefully before use the kit.
- Do not use after expiration date.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use product if the protective external box or the protective aluminium pouches are opened or damaged upon arrival.
- Do not use the test (and positive/negative controls if the final reference contains them) if desiccant material is not present or broken inside the aluminium pouch.
- Components provided in the kit are approved for use in the **CerTest C.*difficile* GDH+Toxin A+B**. Do not mix or use with components from other product kit or lots. Do not use with other commercial kit components (from another supplier).
- All the components of the kits are for single use. Do not reuse.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves or other personal protective equipment such as goggles and mask that will be considered necessary. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous, as well as all the materials that have been exposed to the samples and they must be handled according to the local or national safety regulations. Must be handled in the same manner as an infectious agent. Use proper infection control practices. These practices should include, but are not limited to, personal protective equipment (PPE), such as laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Take necessary precautions during the collection, transport, storage, handling, and disposal of the samples. Each sample must be correctly and unequivocally identified in order to guarantee the correct traceability of the samples.
- A new test must be used for each sample to avoid contaminations errors.
- Clean up spills thoroughly using an appropriate disinfectant.

- Reagents contain preservatives (<0.1% sodium azide). Avoid any contact with skin or mucous membrane. In accordance with Regulation (EC) No 1907/2006 (REACH), **CerTest C.difficile GDH+Toxin A+B** do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration. Material Safety Data Sheet are not included with this device, nevertheless it could be consulted upon requested to CerTest Biotec.
- The certificate of analysis is not included with the device, however, it could be requested to CerTest Biotec in case of need.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Visual interpretation of the results by coloured lines, interpretation of the results must be done by professional user without problems in colour interpretation.
- The tests (vials, devices, positive/negative control if any) should be discarded in a proper biohazard container after testing. These containers should be discarded in accordance with local or national laws or regulations.

7. Specimen collection and preparation


Samples should be collected in a clean container and processed as soon as possible after collection. Follow proper infection control practices.

Samples can be preserved at room temperature if they are going to be analysed into the same day after taking them. The recommendation is processed the sample immediately after collection.

If this is not possible, stool samples could be stored in the refrigerator (2-8°C) for a maximum of 24 hours. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the **stool collection tube** (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 50 mg) and add it to the **stool collection tube**. For liquid samples, add approx. 125 µL in the **stool collection tube** using a micropipette.
2. Close the **stool collection tube**. Shake the tube in order to assure good sample dispersion (3).

Continue immediately with the test procedure. See Annexes. Each sample vial (containing the sample) and cassette with which will be tested that sample should be perfectly recognizable including patient's identification.

8. Limitations of the test

- 1) The test must be carried out within 2 hours after opening the sealed bag.
- 2) An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- 3) The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- 4) **CerTest C. difficile GDH+Toxin A+B** should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- 5) Positive results determine the presence of GDH, Toxin A and Toxin B of *Clostridium difficile* in faecal samples. A positive result should be followed up with additional laboratory techniques (toxigenic culture) to determine the strain. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- 6) A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation persist, a *Clostridium difficile* determination should be carried out, on a sample from an enrichment culture.
- 7) Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool samples whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
- 8) Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

9. Expected values

Clostridium difficile is associated with 95-100% cases of pseudomembranous colitis, 60-75% cases of antibiotic-associated colitis and 35% of cases of antibiotic-associated diarrhoea cases. In addition, nearly 223.900

people in the United States required hospital care for *C. difficile* and at least 12,800 people died in 2017. Although U.S. outpatient antibiotic prescription rate decreased from 2011 to 2016, at least 30% of outpatient antibiotic prescription are estimated to be unnecessary, which highlights the need to improve outpatient prescribing (4–6).

While the development of new treatments is one of these key actions, such investments must be coupled with dedicated efforts toward preventing infections in the first place, slowing the development of resistance through better antibiotic use, and stopping the spread of resistance (6).

10. Performance characteristics

A. Analytical sensitivity (detection limit)

Detection limit value (LoD) of **CerTest C. difficile GDH+Toxin A+B** is 0.39 ng/mL of *Clostridium difficile* antigen GDH, 2 ng/mL of *Clostridium difficile* Toxin A and 3.12 ng/mL of *Clostridium difficile* Toxin B.

There is not any external standard approved as a general and reference standard. The internal standard used for validation of this CerTest product are purified *Clostridium difficile* Glutamate dehydrogenase recombinant protein and purified *Clostridium difficile* Toxin A and Toxin B.

B. Clinical sensitivity and specificity

An external retrospective evaluation (Hospital Clinico Lozano Blesa, Zaragoza) with frozen stool samples was performed using **CerTest C. difficile GDH+Toxin A+B** versus another immunochromatographic (IC) test. Discordant results between the IC tests were analyzed by qPCR. This is the evaluation criteria. The results were as follows:

For Strip A: *Clostridium difficile* antigen GDH

		Evaluation criteria (GDH)		
IC test: CerTest C. difficile GDH+Toxin A+B (GDH)		+	-	Total
	+	86	1	87
	-	3	160	163
	Total	89	161	250

Table 1 CerTest C. difficile GDH+Toxin A+B (GDH) performance versus evaluation criteria.

CerTest C. difficile GDH+Toxin A+B (GDH) vs Evaluation criteria		
	Mean Value	95% confidence interval
Sensitivity	96.6%	90.5-99.3%
Specificity	99.4%	96.6-100.0%
PPV	98.9%	93.8-100.0%
NPV	98.2%	94.7-99.6%

Table 2 CerTest C. difficile GDH+Toxin A+B (GDH) clinical results.

For Strip B: *Clostridium difficile* Toxin A

		Evaluation criteria (Toxin A)		
IC test: CerTest <i>C. difficile</i> GDH+Toxin A+B (Toxin A)		+	-	Total
	+	51	0	51
	-	1	198	199
	Total	52	198	250

Table 3 CerTest *C. difficile* GDH+Toxin A+B (Toxin A) performance versus evaluation criteria.

CerTest <i>C. difficile</i> GDH+Toxin A+B (Toxin A) vs Evaluation criteria		
	Mean Value	95% confidence interval
Sensitivity	98.1%	89.7-100.0%
Specificity	100.0%	98.2-100.0%
PPV	100.0%	93.0-100.0%
NPV	99.5%	97.2-100.0%

Table 4 CerTest *C. difficile* GDH+Toxin A+B (Toxin A) clinical results.**For Strip C: *Clostridium difficile* Toxin B**

		Evaluation criteria (Toxin B)		
IC test: CerTest <i>C. difficile</i> GDH+Toxin A+B (Toxin B)		+	-	Total
	+	40	0	40
	-	4	206	210
	Total	44	206	250

Table 5 CerTest *C. difficile* GDH+Toxin A+B (Toxin B) performance versus evaluation criteria.

CerTest <i>C. difficile</i> GDH+Toxin A+B (Toxin B) vs Evaluation criteria		
	Mean Value	95% confidence interval
Sensitivity	90.9%	78.3-97.5%
Specificity	100.0%	98.2-100.0%
PPV	100.0%	91.2-100.0%
NPV	98.1%	95.2-99.5%

Table 6 CerTest *C. difficile* GDH+Toxin A+B (Toxin B) clinical results.

The most recent external retrospective evaluation performed to the test (with frozen stool samples) showed the following results:

- ***Clostridium difficile* antigen GDH:** sensitivity 96.6% (95% confidence interval: 90.5-99.3%) and specificity 99.4% (95% confidence interval: 96.6-100.0%).
- **Toxin A:** Sensitivity 98.1%, (95% confidence interval: 89.7-100.0%) and specificity 100.0% (95% confidence interval: 98.2-100.0%).

- **Toxin B:** Sensitivity 90.9%, (95% confidence interval: 78.3-97.5%) and specificity 100.0% (95% confidence interval: 98.2-100.0%).

The results showed a high sensitivity and specificity to detect *Clostridium difficile* antigen GDH, Toxin A and/or Toxin B of *Clostridium difficile* using **CerTest C. difficile GDH+Toxin A+B**.

For Strip A: *Clostridium difficile* antigen GDH

Nº	Site	Sample type	Reference	Target
1	Hospital Clínico Lozano Blesa (Zaragoza, Spain) April-May 2021	Human stool samples	Code (17.1.1)	<i>Clostridium difficile</i> antigen GDH
2	Hospital Clínico Lozano Blesa (Zaragoza, Spain) April-June 2017	Human stool samples	Code (17.1.2)	<i>Clostridium difficile</i> antigen GDH
3	Louvain, (France) January- December 2016	Human stool samples	Code (17.1.3)	<i>Clostridium difficile</i> antigen GDH
4	ECCMID April 2016, Amsterdam	Human stool samples	Code (17.1.4)	<i>Clostridium difficile</i> antigen GDH

Table 7 Resume of Clinical evaluations. Site, sample type, target.

Nº	TP	FP	FN	TN	N	Sens. (%) CI (%)	Spec. (%) CI (%)	PPV (%) CI (%)	NPV (%) CI (%)
1	86	1	3	160	250	96.6 (90.5-99.3)	99.4 (96.6-100.0)	98.9 (93.8-100.0)	98.2 (94.7-99.6)
2	51	0	0	75	126	100.0 (93.0-100.0)	100.0 (95.2-100.0)	100.0 (93.0-100.0)	100.0 (95.2-100.0)
3	130	8	8	332	478	94.2 (88.9-97.5)	97.6 (95.4-99.0)	94.2 (88.9-97.5)	97.6(95.4-99.0)
4	27	8	6	117	158	81.8 (64.5-93.0)	93.6 (87.8-97.2)	77.1 (59.9-89.6)	95.1 (89.7-98.2)
Total	294	17	17	684	1012	94.5 (91.4-96.8)	97.6 (96.1-98.6)	94.5 (91.4-96.8)	97.6 (96.1-98.6)

Table 8 Clinical evaluations and global results. True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Predictive Positive Values (PPV), Predictive Negative Values (NPV), 95% interval of confidence (CI) for CerTest C. *difficile* GDH+Toxin A+B (GDH).

For Strip B: *Clostridium difficile* Toxin A

Ref.	Location	Sample type	Reference	Target
1	Hospital Clínico Lozano Blesa (Zaragoza, Spain) April-June 2017	Human stool samples	Code (17.1.1)	<i>Clostridium difficile</i> Toxin A
2	Hospital Clínico Lozano Blesa (Zaragoza, Spain) April-May 2021	Human stool samples	Code (17.1.2)	<i>Clostridium difficile</i> Toxin A

Table 9 Resume of Clinical evaluations. Location, sample type, target.

Ref.	TP	FP	FN	TN	N	Sens. (%) CI (%)	Spec. (%) CI (%)	PPV (%) CI (%)	NPV (%) CI (%)
1	28	0	1	97	126	96.9 (92.2-99.9)	100.0 (96.2-100.0)	100.0 (87.7-100.0)	99.0 (94.4-100.0)
2	51	0	1	198	250	98.1 (89.7-100.0)	100.0 (98.2-100.0)	100.0 (93.0-100.0)	99.5 (97.2-100.0)
Total	79	0	2	295	376	97.5 (91.4-99.7)	100.0 (98.8-100.0)	100.0 (95.4-100.0)	99.3 (97.6-99.9)

Table 10 Clinical evaluations and global results. True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Predictive Positive Values (PPV), Predictive Negative Values (NPV), 95% interval of confidence (CI) for CerTest *C. difficile* GDH+Toxin A+B (Toxin A).

For Strip C: *Clostridium difficile* Toxin B

Ref.	Location	Sample type	Reference	Target
1	Hospital Clínico Lozano Blesa (Zaragoza, Spain) April-June 2017	Human stool samples	Code (17.1.1.)	<i>Clostridium difficile</i> Toxin B
2	Hospital Clínico Lozano Blesa (Zaragoza, Spain) April-May 2021	Human stool samples	Code (17.1.2)	<i>Clostridium difficile</i> Toxin B

Table 11 Resume of Clinical evaluations. Location, sample type, target.

Ref.	TP	FP	FN	TN	N	Sens. (%) CI (95%)	Spec. (%) CI (%)	PPV (%) CI (%)	NPV (%) CI (%)
1	36	1	0	89	126	100.0 (90.3-100.0)	98.9 (94.0-100.0)	97.3 (85.8-99.9)	100.0 (95.9-100.0)
2	40	0	4	206	250	90.9 (78.3-97.5)	100.0 (98.2-100.0)	100.0 (91.2-100.0)	98.1 (95.2-99.5)
Total	76	1	4	295	376	95.0 (87.7-98.6)	99.7 (98.1-100.0)	98.7 (93.0-100.0)	98.7 (96.6-99.6)

Table 12 Clinical evaluations and global results. True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Predictive Positive Values (PPV), Predictive Negative Values (NPV), 95% interval of confidence (CI) for CerTest *C. difficile* GDH+Toxin A+B (Toxin B).

C. Hook Effect

CerTest C. difficile GDH+Toxin A+B does not show inhibitory hook effect at concentration higher than $10^5 \times \text{LoD}$ for *Clostridium difficile* antigen GDH, higher than $10^2 \times \text{LoD}$ for *C. difficile* Toxin A and higher than $10^4 \times \text{LoD}$ for *C. difficile* Toxin B.

D. Cross Reactivity

An evaluation was performed to determine the cross reactivity of **CerTest C. difficile GDH+Toxin A+B**, no cross reactivity against organism, pathogens, substances that could cause infections:

For Strip A: *Clostridium difficile* antigen GDH

Cross reactivity									
Virus	Adenovirus	-	Coronavirus	-	Influenza B	-	Norovirus GII	-	Rotavirus
	Astrovirus	-	Influenza A	-	Norovirus GI	-	Respiratory Syncytial Virus	-	
Bacteria	<i>Campylobacter coli</i>	-	<i>Clostridium Haemolyticum</i>	-	<i>Helicobacter pylori</i>	-	<i>Salmonella typhimurium</i>	-	<i>Streptococcus pyogenes</i>
	<i>Campylobacter jejuni</i>	-	<i>Clostridium Novyi</i>	-	<i>Legionella pneumophila</i>	-	<i>Shigella boydii</i>	-	<i>Yersinia enterocolitica</i> O:3
	<i>Clostridium difficile</i> Toxin A	-	<i>Clostridium Tetani</i>	-	<i>Listeria monocytogenes</i>	-	<i>Shigella dysenteriae</i>	-	<i>Yersinia enterocolitica</i> O:9
	<i>Clostridium difficile</i> Toxin B	-	<i>Clostridium Septicum</i>	-	<i>Peptostreptococcus anaerobius</i>	-	<i>Shigella flexneri</i>	-	
	<i>Clostridium perfringens</i>	-	<i>Escherichia coli</i> O:157	-	<i>Salmonella enteritidis</i>	-	<i>Shigella sonnei</i>	-	
	<i>Clostridium bifermentans</i>	-	<i>Escherichia coli</i> O:026	-	<i>Salmonella paratyphi</i>	-	<i>Staphylococcus aureus</i>	-	
	<i>Clostridium Butyricum</i>	-	<i>Escherichia coli</i> O:111	-	<i>Salmonella typhi</i>	-	<i>Streptococcus pneumococcal</i>	-	
Parasites	<i>Cryptosporidium Parvum</i>	-	<i>Entamoeba histolytica</i>	-	<i>Giardia lamblia</i>	-			
Others	Pig haemoglobin	-	Bovine haemoglobin	-	Bovine Lactoferrin	-	Bovine Transferrin	-	

Table 13 Cross reactivity results (GDH).

A specificity assay was performed for **CerTest *Clostridium difficile* antigen GDH**. This concluded that the test, could detect the following antigens: *C. sporogenes* (CECT 485) and *C. botulinum* (CECT 551).

In this specificity assay, *Clostridium sordellii* (ATCC 9714) was evaluated and no cross reactivity was found. In few intercomparison studies **CerTest *Clostridium difficile* antigen GDH** have a positive signal with this pathogen, but the nature of this positive result could not be established.

For Strip B: *Clostridium difficile* Toxin A

Cross reactivity									
Virus	Adenovirus	-	Coronavirus	-	Influenza B	-	Norovirus GII	-	Rotavirus
	Astrovirus	-	Influenza A	-	Norovirus GI	-	Respiratory Syncytial Virus	-	
Bacteria	<i>Campylobacter coli</i>	-	<i>Escherichia coli</i> O:157	-	<i>Listeria monocytogenes</i>	-	<i>Shigella boydii</i>	-	<i>Streptococcus pneumococcal</i>
	<i>Campylobacter jejuni</i>	-	<i>Escherichia coli</i> O:026	-	<i>Salmonella enteritidis</i>	-	<i>Shigella dysenteriae</i>	-	<i>Streptococcus pyogenes</i>
	<i>Clostridium difficile</i> antigen GDH	-	<i>Escherichia coli</i> O:111	-	<i>Salmonella paratyphi</i>	-	<i>Shigella flexneri</i>	-	<i>Yersinia enterocolitica</i> O:3
	<i>Clostridium difficile</i> Toxin B	-	<i>Helicobacter pylori</i>	-	<i>Salmonella typhi</i>	-	<i>Shigella sonnei</i>	-	<i>Yersinia enterocolitica</i> O:9
	<i>Clostridium perfringens</i>	-	<i>Legionella pneumophila</i>	-	<i>Salmonella typhimurium</i>	-	<i>Staphylococcus aureus</i>	-	
Parasites	<i>Cryptosporidium Parvum</i>	-	<i>Entamoeba histolytica</i>	-	<i>Entamoeba dispar</i>	-	<i>Giardia lamblia</i>	-	
Others	Pig haemoglobin	-	Bovine haemoglobin	-					

Table 14 Cross reactivity results (Toxin A).

For Strip C: *Clostridium difficile* Toxin B

Cross reactivity									
Virus	Adenovirus	-	Coronavirus	-	Influenza B	-	Norovirus GII	-	Rotavirus
	Astrovirus	-	Influenza A	-	Norovirus GI	-	Respiratory Syncytial Virus	-	
Bacteria	<i>Campylobacter coli</i>	-	<i>Escherichia coli</i> O:157	-	<i>Listeria monocytogenes</i>	-	<i>Shigella boydii</i>	-	<i>Streptococcus pneumococcal</i>
	<i>Campylobacter jejuni</i>	-	<i>Escherichia coli</i> O:026	-	<i>Salmonella enteritidis</i>	-	<i>Shigella dysenteriae</i>	-	<i>Streptococcus pyogenes</i>
	<i>Clostridium difficile</i> antigen GDH	-	<i>Escherichia coli</i> O:111	-	<i>Salmonella paratyphi</i>	-	<i>Shigella flexneri</i>	-	<i>Yersinia enterocolitica</i> O:3
	<i>Clostridium difficile</i> Toxin A	-	<i>Helicobacter pylori</i>	-	<i>Salmonella typhi</i>	-	<i>Shigella sonnei</i>	-	<i>Yersinia enterocolitica</i> O:9
	<i>Clostridium perfringens</i>	-	<i>Legionella pneumophila</i>	-	<i>Salmonella typhimurium</i>	-	<i>Staphylococcus aureus</i>	-	
Parasites	<i>Cryptosporidium Parvum</i>	-	<i>Entamoeba histolytica</i>	-	<i>Entamoeba dispar</i>	-	<i>Giardia lamblia</i>	-	
Others	Pig haemoglobin	-	Bovine haemoglobin	-					

Table 15 Cross reactivity results (Toxin B).

E. Interferences

An evaluation was performed to determine the possible interferences of **CerTest C. difficile GDH+Toxin A+B**; no interferences against the substances tested were detected:

For Strip A: *Clostridium difficile* antigen GDH

	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
Exogenous	Metronidazole	3.0 mg/mL	No interference	Loperamide hydrochloride (Fortasec)	0.15 mg/mL	No interference
	Ampicillin	3.0 mg/mL	No interference	Heparin (Hibor)	350.0 IU/mL	No interference
	Oseltamivir	0.3 mg/mL	No interference	Almagato (Almax)	3.0 mg/mL	No interference
	Amantadine	0.3 mg/mL	No interference	Fosfamycin (Monurol)	3.0·10 ⁻³ mg/mL	No interference
	Ribavirin	3.0 mg/mL	No interference	Acetylcysteine (Fluimucil)	3.0 mg/mL	No interference
	Codeine (Toseina)	0.2 mg/mL	No interference	Dexketoprofen trometamol (Enantyum)	0.3 mg/mL	No interference
	Benzocaine (Angileptol)	0.3 mg/mL	No interference	Levofloxacin	3.0 mg/mL	No interference
	Cloperastine (Flutox)	0.3 mg/mL	No interference	Ciprofloxacin	0.3 mg/mL	No interference
	Carbocisteine (Iniston mucolítico)	3.0 mg/mL	No interference	Rifampicin (Rifaldin)	0.3 mg/mL	No interference
	Loratadine	0.3 mg/mL	No interference	Phenoxymethylpenicillin potassium	3.0 mg/mL	No interference
	Dexchlorpheniramine (Polaramine)	0.3 mg/mL	No interference	Ambroxol hydrochloride (Mucosan)	0.3 mg/mL	No interference
	Ebastine (Ebastel)	3.0 mg/mL	No interference	Macrogol 3350 (Movicol)	3.0 mg/mL	No interference
	Acetyl Salicylic (Adiro)	0.3 mg/mL	No interference	Lysine Carbocysteinate (Pectox)	3.0·10 ⁻² mg/mL	No interference
	Ibuprofen (Espidifen)	3.0 mg/mL	No interference	Hydroxyzine dihydrochloride	0.3 mg/mL	No interference
	Paracetamol (Dolocatil)	5.0 mg/mL	No interference	Lorazepam	3.0·10 ⁻³ mg/mL	No interference
	Metamizole (Nolotil)	5.0 mg/mL	No interference	Amoxicillin	3.0 mg/mL	No interference
	Prednisone	0.3 mg/mL	No interference	Mercaptopurine	0.3 mg/mL	No interference
	Omeprazole	2.0·10 ⁻³ mg/mL	No interference	Biotine	0.10 mg/mL	No interference
	Naso GEL	0.9 mg/mL	No interference	Sore Throat Phenol spray	0.5 mg/mL	No interference
	CVS Nasal Spray (Cromolyn)	4 mg/mL	No interference	Tobramycin	0.3 mg/mL	No interference
	Afrin (Oxymetazoline)	5.0·10 ⁻² mg/mL	No interference	Mupirocin	2.5·10 ⁻² mg/mL	No interference
	CVS Nasal Drops (Phenylephrine)	10.0 mg/mL	No interference	Fluticasone Propionate	5.0·10 ⁻² mg/mL	No interference
	ZICAM	1.0 mg/mL	No interference	Homeopathic	DIL 1/10	No interference
	Phenylpropanolamine	0.3 mg/mL	No interference			
Endogenous	Human haemoglobin	5.0 µg/mL	No interference	Human lactoferrin	5.0 µg/mL	No interference
	Human transferrin	0.5 µg/mL	No interference	Mucine	0.5 %	No interference
	Human calprotectin	5.0 µg/mL	No interference	Human blood	5.0 %	No interference

Table 16 *Clostridium difficile* antigen GDH interferences results.

For Strip B: *Clostridium difficile* Toxin A

	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
Exogenous	Metronidazole	3.0 mg/mL	No interference	Loperamide hydrochloride (<i>Fortasec</i>)	0.15 mg/mL	No interference
	Ampicillin	3.0 mg/mL	No interference	Heparin (<i>Hibor</i>)	350.0 IU/mL	No interference
	Oseltamivir	0.3 mg/mL	No interference	Almagato (<i>Almax</i>)	3.0 mg/mL	No interference
	Amantadine	0.3 mg/mL	No interference	Fosfomicin (<i>Monurol</i>)	3.0·10 ⁻³ mg/mL	No interference
	Ribavirin	3.0 mg/mL	No interference	Acetylcysteine (<i>Fluimucil</i>)	3.0 mg/mL	No interference
	Codeine (<i>Toseina</i>)	0.2 mg/mL	No interference	Dexketoprofen trometamol (<i>Enantyum</i>)	0.3 mg/mL	No interference
	Benzocaine (<i>Angileptol</i>)	0.3 mg/mL	No interference	Levofloxacin	3.0 mg/mL	No interference
	Cloperastine (<i>Flutox</i>)	0.3 mg/mL	No interference	Ciprofloxacin	0.3 mg/mL	No interference
	Carbocisteine (<i>Iniston mucolítico</i>)	3.0 mg/mL	No interference	Rifampicin (<i>Rifaldin</i>)	0.3 mg/mL	No interference
	Loratadine	0.3 mg/mL	No interference	Phenoxyethylpenicillin potassium	3.0 mg/mL	No interference
	Dexchlorpheniramine (<i>Polaramine</i>)	0.3 mg/mL	No interference	Ambroxol hydrochloride (<i>Mucosan</i>)	0.3 mg/mL	No interference
	Ebastine (<i>Ebastel</i>)	3.0 mg/mL	No interference	Macrogol 3350 (<i>Movicol</i>)	3.0 mg/mL	No interference
	Acetyl Salicylic (<i>Adiro</i>)	0.3 mg/mL	No interference	Lysine Carbocysteinate (<i>Pectox</i>)	3.0·10 ⁻² mg/mL	No interference
	Ibuprofen (<i>Espidifen</i>)	3.0 mg/mL	No interference	Hydroxyzine dihydrochloride	0.3 mg/mL	No interference
	Paracetamol (<i>Dolocatil</i>)	5.0 mg/mL	No interference	Lorazepam	3.0·10 ⁻³ mg/mL	No interference
	Metamizole (<i>Nolotil</i>)	5.0 mg/mL	No interference	Amoxicillin	3.0 mg/mL	No interference
	Prednisone	0.3 mg/mL	No interference	Mercaptopurine	0.3 mg/mL	No interference
	Omeprazole	2.0·10 ⁻³ mg/mL	No interference	Biotine	0.10 mg/mL	No interference
	Naso GEL	0.9 mg/mL	No interference	Sore Throat Phenol spray	0.5 mg/mL	No interference
	CVS Nasal Spray (<i>Cromolyn</i>)	4.0 mg/mL	No interference	Tobramycin	0.3 mg/mL	No interference
	Afrin (<i>Oxymetazoline</i>)	5.0·10 ⁻² mg/mL	No interference	Mupirocin	2.5·10 ⁻² mg/mL	No interference
	CVS Nasal Drops (<i>Phenylephrine</i>)	10.0 mg/mL	No interference	Fluticasone Propionate	5.0·10 ⁻² mg/mL	No interference
	ZICAM	1.0 mg/mL	No interference	Homeopathic	DIL 1/10	No interference
	Phenylpropanolamine	0.3 mg/mL	No interference			
Endogenous	Human haemoglobin	5.0 µg/mL	No interference	Human lactoferrin	5.0 µg/mL	No interference
	Human transferrin	0.5 µg/mL	No interference	Mucine	0.5 %	No interference
	Human calprotectin	5.0 µg/mL	No interference	Human blood	5.0 %	No interference

Table 11 *Clostridium difficile* toxin A interferences results.

For Strip C: *Clostridium difficile* Toxin B

	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
Exogenous	Metronidazole	3.0 mg/mL	No interference	Loperamide hydrochloride (<i>Fortasec</i>)	0.15 mg/mL	No interference
	Ampicillin	3.0 mg/mL	No interference	Heparin (<i>Hibor</i>)	350.0 IU/mL	No interference
	Oseltamivir	0.3 mg/mL	No interference	Almagato (<i>Almax</i>)	3.0 mg/mL	No interference
	Amantadine	0.3 mg/mL	No interference	Fosfomicin (<i>Monurol</i>)	3.0·10 ⁻³ mg/mL	No interference
	Ribavirin	3.0 mg/mL	No interference	Acetylcysteine (<i>Fluimucil</i>)	3.0 mg/mL	No interference
	Codeine (<i>Toseina</i>)	0.2 mg/mL	No interference	Dexketoprofen trometamol (<i>Enantyum</i>)	0.3 mg/mL	No interference
	Benzocaine (<i>Angileptol</i>)	0.3 mg/mL	No interference	Levofloxacin	3.0 mg/mL	No interference
	Cloperastine (<i>Flutox</i>)	0.3 mg/mL	No interference	Ciprofloxacin	0.3 mg/mL	No interference
	Carbocisteine (<i>Iniston mucolítico</i>)	3.0 mg/mL	No interference	Rifampicin (<i>Rifaldin</i>)	0.3 mg/mL	No interference
	Loratadine	0.3 mg/mL	No interference	Phenoxymethylpenicillin potassium	3.0 mg/mL	No interference
	Dexchlorpheniramine (<i>Polaramine</i>)	0.3 mg/mL	No interference	Ambroxol hydrochloride (<i>Mucosan</i>)	0.3 mg/mL	No interference
	Ebastine (<i>Ebastel</i>)	3.0 mg/mL	No interference	Macrogol 3350 (<i>Movicol</i>)	3.0 mg/mL	No interference
	Acetyl Salicylic (<i>Adiro</i>)	0.3 mg/mL	No interference	Lysine Carbocysteinate (<i>Pectox</i>)	3.0·10 ⁻² mg/mL	No interference
	Ibuprofen (<i>Espidifen</i>)	3.0 mg/mL	No interference	Hydroxyzine dihydrochloride	0.3 mg/mL	No interference
	Paracetamol (<i>Dolocatil</i>)	5.0 mg/mL	No interference	Lorazepam	3.0·10 ⁻³ mg/mL	No interference
	Metamizole (<i>Nolotil</i>)	5.0 mg/mL	No interference	Amoxicillin	3.0 mg/mL	No interference
	Prednisone	0.3 mg/mL	No interference	Mercaptopurine	0.3 mg/mL	No interference
	Omeprazole	2.0·10 ⁻³ mg/mL	No interference	Biotine	0.10 mg/mL	No interference
	Naso GEL	0.9 mg/mL	No interference	Sore Throat Phenol spray	0.5 mg/mL	No interference
	CVS Nasal Spray (<i>Cromolyn</i>)	4.0 mg/mL	No interference	Tobramycin	0.3 mg/mL	No interference
	Afrin (<i>Oxymetazoline</i>)	5.0·10 ⁻² mg/mL	No interference	Mupirocin	2.5·10 ⁻² mg/mL	No interference
	CVS Nasal Drops (<i>Phenylephrine</i>)	10.0 mg/mL	No interference	Fluticasone Propionate	5.0·10 ⁻² mg/mL	No interference
	ZICAM	1.0 mg/mL	No interference	Homeopathic	DIL 1/10	No interference
	Phenylpropanolamine	0.3 mg/mL	No interference			
Endogenous	Human haemoglobin	5.0 µg/mL	No interference	Human lactoferrin	5.0 µg/mL	No interference
	Human transferrin	0.5 µg/mL	No interference	Mucine	0.5 %	No interference
	Human calprotectin	5.0 µg/mL	No interference	Human blood	5.0 %	No interference

Table 12 *Clostridium difficile* toxin B interferences results.**F. Repeatability and Reproducibility**

Following the study of repeatability and reproducibility of **CerTest C. *difficile* GDH+Toxin A+B**, performed using different internal samples, negative and positive, no significative differences have been observed within the evaluations.

ANNEX 1

KIT 20 CARD C. DIFFICILE GDH+TOXIN A+B. Catalogue reference: GX872001V

Annex for the following reference:

PRODUCT	REFERENCE
CerTest C. difficile GDH+Toxin A+B Combo Card (20 tests kit with stool vials)	GX872001V

Table A 1. Reference for kit 20 combo card of CerTest C.difficile GDH+Toxin A+B.

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
- 20x CerTest C. difficile GDH+Toxin A+B (combo card test format) - 20x Stool collection tubes with diluent - 1x IFU Instructions for use (including Quick guide use).	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Specimen collection container


Table A 1.1 Materials necessities for test performance.

II. Test procedures (samples)

Follow proper infection control practices. Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the **stool collection tube** in order to assure good sample dispersion.
2. Remove **CerTest C. difficile GDH+Toxin A+B** (combo card test) from its sealed bag just before using it.
3. Take the **stool collection tube**, cut the end of the cap (1) and dispense 3 drops in the circular window marked with the letter A (2), 3 drops, using the same tube, in the circular window marked with letter B (3) and 3 drops, using the same tube, in the circular window marked with letter C (4). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.





If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

III. Interpretation of results

(please, refer to the illustration below)

	GDH	Toxin A	Toxin B		Interpretation of the results
1.	NEGATIVE	NEGATIVE	NEGATIVE		There are no GDH, Toxin A or Toxin B presence. No infection caused by <i>Clostridium difficile</i>. A negative result should be followed up with additional laboratory techniques to confirm the results.
	GREEN	GREEN	GREEN		
2.	POSITIVE	POSITIVE	POSITIVE		There are GDH, Toxin A and Toxin B presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result.
	GREEN-RED	GREEN-RED	GREEN-RED		

	GDH	Toxin A	Toxin B		Interpretation of the results
					A positive result should be followed up with additional laboratory techniques to confirm the results.
3.	POSITIVE	POSITIVE	NEGATIVE		There is GDH and Toxin A presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN-RED	GREEN		
4.	POSITIVE	NEGATIVE	POSITIVE		There is GDH and Toxin B presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN	GREEN-RED		
5.	POSITIVE	NEGATIVE	NEGATIVE		There is GDH presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN	GREEN		
6.	NEGATIVE	POSITIVE	POSITIVE		If this result appears it must be repeat the test using a fresh sample. If results are again positive for Toxin A and B and negative for GDH, the sample should be considered positive for Toxin A and B. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result.
	GREEN	GREEN-RED	GREEN-RED		



	GDH	Toxin A	Toxin B		Interpretation of the results
					A positive result should be followed up with additional laboratory techniques to confirm the results.
7.	NEGATIVE	POSITIVE	NEGATIVE		If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin A and negative for GDH, the sample should be considered positive for Toxin A. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
		GREEN	GREEN-RED		
8.	NEGATIVE	NEGATIVE	POSITIVE		If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin B and negative for GDH, the sample should be considered positive for Toxin B. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
		GREEN	GREEN		
9.	INVALID				Invalid result either A, B or C, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid.
	Any other result				

Table A 1.2. CerTest C. difficile GDH+Toxin A+B performance interpretation of the results.

Whatever results detailed between 6 and 8, are not very common. A very low percentage of specimens might result negative for GDH but positive for toxins.

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality control

INTERNAL CONTROL: control line

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

ANNEX 2

KIT 20 CARD C. *DIFFICILE* GDH+TOXIN A+B. Catalogue reference: GX872001VC

Annex for the following reference:

PRODUCT	REFERENCE
CerTest C. <i>difficile</i> GDH+Toxin A+B Combo Card (20 tests kit with stool vials + PC)	GX872001VC

Table A 2. Reference for kit 20 combo card of CerTest C. *difficile* GDH+Toxin A+B with positive control.

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
- 20x CerTest C. <i>difficile</i> GDH+Toxin A+B (combo card test format) - 20x Stool collection tubes with diluent -1x CerTest C. <i>difficile</i> antigen GDH Positive Control -1x CerTest C. <i>difficile</i> Toxin A Positive Control -1x CerTest C. <i>difficile</i> Toxin B Positive Control -3x Testing tube -3x Disposable pipette -3x Positive control reagent tubes (red cap) -1x IFU Instructions for use (including Quick guide use).	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Specimen collection container


Table A 2.1 Materials necessities for test performance.

II. Test procedures (samples)


Follow proper infection control practices. Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the **stool collection tube** in order to assure good sample dispersion.
2. Remove **CerTest C. difficile GDH+Toxin A+B** (combo card test) from its sealed bag just before using it.
3. Take the **stool collection tube**, cut the end of the cap (1) and dispense 3 drops in the circular window marked with the letter A (2), 3 drops, using the same tube, in the circular window marked with letter B (3) and 3 drops, using the same tube, in the circular window marked with letter C (4). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.


(2) Add 3 drops in the circular window (A)

(3) Add 3 drops in the circular window (B)





(4) Add 3 drops in the circular window (C)





Wait 10 minutes

III. Interpretation of results

(please, refer to the illustration below)

	GDH	Toxin A	Toxin B		Interpretation of the results
1.	NEGATIVE	NEGATIVE	NEGATIVE		There are no GDH, Toxin A or Toxin B presence. No infection caused by <i>Clostridium difficile</i>. A negative result should be followed up with additional laboratory techniques to confirm the results.
	GREEN	GREEN	GREEN		
2.	POSITIVE	POSITIVE	POSITIVE		There are GDH, Toxin A and Toxin B presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN-RED	GREEN-RED		
3.	POSITIVE	POSITIVE	NEGATIVE		There is GDH and Toxin A presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN-RED	GREEN		
4.	POSITIVE	NEGATIVE	POSITIVE		There is GDH and Toxin B presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN	GREEN-RED		

	GDH	Toxin A	Toxin B		Interpretation of the results
5.	POSITIVE	NEGATIVE	NEGATIVE		There is GDH presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN	GREEN		
6.	NEGATIVE	POSITIVE	POSITIVE		If this result appears it must be repeat the test using a fresh sample. If results are again positive for Toxin A and B and negative for GDH, the sample should be considered positive for Toxin A and B. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN	GREEN-RED	GREEN-RED		
7.	NEGATIVE	POSITIVE	NEGATIVE		If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin A and negative for GDH, the sample should be considered positive for Toxin A. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN	GREEN-RED	GREEN		
8.	NEGATIVE	NEGATIVE	POSITIVE		If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin B and negative for GDH, the sample should be considered positive for Toxin B. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line),
	GREEN	GREEN	GREEN-RED		

	GDH	Toxin A	Toxin B		Interpretation of the results
					independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
9.	INVALID				Invalid result either A, B or C, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid.
	Any other result				

Table A 2.2 CerTest *C. difficile* GDH+Toxin A+B performance interpretation of the results.

Whatever results detailed between 6 and 8, are not very common. A very low percentage of specimens might result negative for GDH but positive for toxins.

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality control

CerTest *C. difficile* GDH+Toxin A+B (combo card test) presents four controls:

INTERNAL CONTROL: control line

Internal control: internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

EXTERNAL CONTROL: POSITIVE CONTROL SWABS

External control: CerTest C. difficile antigen GDH Positive Control, CerTest C. difficile Toxin A Positive Control and CerTest C. difficile Toxin B Positive Control (positive control swabs).

A. External positive controls

It is recommended to use the same device for both external positive controls.

1. CerTest C. difficile antigen GDH Positive Control


CerTest C. difficile antigen GDH Positive Control is an external quality control for **CerTest C. difficile GDH+Toxin A+B** (Strip A). Use **CerTest C. difficile antigen GDH Positive Control** to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C, **CerTest C. difficile antigen GDH Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest C. difficile antigen GDH Positive Control** must remain in the sealed pouch until use. Do not freeze.

1.1. Positive control procedure

Allow tests and CerTest C. difficile antigen GDH Positive Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Take the **positive control reagent tube**, cut the end of the cap.
2. Add 15 drops (1) from the **positive control reagent tube** into the **testing tube** and immediately introduce **CerTest C. difficile antigen GDH Positive Control**.
3. Mix the solution by rotating **CerTest C. difficile antigen GDH Positive Control** forcefully against the side of the **testing tube** at least 15 seconds. Extract as much liquid as possible from the **CerTest C. difficile antigen GDH Positive Control**, squeezing the sides of the **testing tube** as **CerTest C. difficile antigen GDH Positive Control** is withdrawn (2). Discard **CerTest C. difficile antigen GDH Positive Control**.
4. Remove the **CerTest C. difficile GDH+Toxin A+B** (combo card format) from its sealed bag just before using it.
5. Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** in the circular window marked with the letter A (3).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

2. CerTest C. difficile Toxin A Positive Control

CerTest C. difficile Toxin A Positive Control is an external quality control for **CerTest C. difficile GDH+Toxin A+B** (Strip B). Use **CerTest C. difficile Toxin A Positive Control** to check that the extraction reagents and the test are working properly.


Store as packaged in the sealed pouch at 2-30°C, **CerTest C. difficile Toxin A Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest C. difficile Toxin A Positive Control** must remain in the sealed pouch until use. Do not freeze.

2.1. Positive control procedure

Allow tests and CerTest C. difficile Toxin A Positive Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Take the **positive control reagent tube**, cut the end of the cap.
2. Add 15 drops (1) from **positive control reagent tube** into the **testing tube** and immediately introduce **CerTest C. difficile Toxin A Positive Control**.
3. Mix the solution by rotating **CerTest C. difficile Toxin A Positive Control** forcefully against the side of the **testing tube** at least 15 seconds. Extract as much liquid as possible from **CerTest C. difficile Toxin A Positive Control**, squeezing the sides of the **testing tube** as **CerTest C. difficile Toxin A Positive Control** is withdrawn (2). Discard **CerTest C. difficile Toxin A Positive Control**.
4. Remove the **CerTest C. difficile GDH+Toxin A+B** (combo card format) from its sealed bag just before using it.
5. Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** in the circular window marked with the letter B (3).

6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

3. CerTest C. difficile Toxin B Positive Control

CerTest C. difficile Toxin B Positive Control is an external quality control for **CerTest C. difficile GDH+Toxin A+B** (Strip C). Use **CerTest C. difficile Toxin B Positive Control** to check that the extraction reagents and the test are working properly.


Store as packaged in the sealed pouch at 2-30°C, **CerTest C. difficile Toxin B Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest C. difficile Toxin B Positive Control** must remain in the sealed pouch until use. Do not freeze.

3.1. Positive control procedure

Allow tests and CerTest C. difficile Toxin B Positive Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Take the **positive control reagent tube**, cut the end of the cap.
2. Add 15 drops (1) from the **positive control reagent tube** into the **testing tube** and immediately introduce the **CerTest C. difficile Toxin B Positive Control**.
3. Mix the solution by rotating the **CerTest C. difficile Toxin B Positive Control** forcefully against the side of the **Testing tube** at least 15 seconds. Extract as much liquid as possible from the **CerTest C. difficile Toxin B Positive Control**, squeezing the sides of the **testing tube** as the **CerTest C. difficile Toxin B Positive Control** is withdrawn (2). Discard **CerTest C. difficile Toxin B Positive Control**.
4. Remove the **CerTest C. difficile GDH+Toxin A+B** (combo card format) from its sealed bag just before using it.
5. Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** in the circular window marked with the letter C (3).



6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Wait 10 minutes

Positive controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

A.1 Interpretation positive control procedure

	CerTest C. <i>difficile</i> antigen GDH Positive Control	CerTest C. <i>difficile</i> Toxin A Positive Control	CerTest C. <i>difficile</i> Toxin B Positive Control		Interpretation of the results
1.	POSITIVE (One device for three positive controls)				A GREEN line (control line (C)) and a RED line (test line (T)) in the three strips, in the results window, appear at the same time due to presence of <i>C.difficile</i> antigens in CerTest C. difficile antigen GDH Positive Control , in strip A, in the CerTest C. difficile Toxin A Positive Control , in strip B, and in CerTest C. difficile Toxin B Positive Control , in strip C.
	GREEN-RED	GREEN-RED	GREEN-RED		
2.	POSITIVE (One device for CerTest C. difficile antigen GDH Positive Control)				A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window in strip A due to presence of antigen GDH in CerTest C. difficile antigen GDH Positive Control .
	GREEN-RED	-	-		



	CerTest C. <i>difficile</i> antigen GDH Positive Control	CerTest C. <i>difficile</i> Toxin A Positive Control	CerTest C. <i>difficile</i> Toxin B Positive Control		Interpretation of the results
3.	POSITIVE (One device for CerTest C. <i>difficile</i> Toxin A Positive Control)				A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window in strip B due to presence of Toxin A in CerTest C. <i>difficile</i> Toxin A Positive Control.
	-	GREEN-RED	-		
4.	POSITIVE (One device for CerTest C. <i>difficile</i> Toxin B Positive Control)				A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window in strip C due to presence of Toxin B in CerTest C. <i>difficile</i> Toxin B Positive Control.
	-	-	GREEN-RED		
4.	INVALID			Invalid result either A, B or C, we recommend repeating the assay using the same sample with another test.	Total absence of the control-coloured line (GREEN), regardless the appearance or not of the test line (RED), or only the control (C) coloured line (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new CerTest C. <i>difficile</i> antigen GDH Positive Control, CerTest C. <i>difficile</i> Toxin A Positive Control and CerTest C. <i>difficile</i> Toxin B Positive Control. If the situation still persists, discontinue using the test kit and contact your local distributor. Any result read after more than 10 minutes will be considered invalid.
	Any other result				

Table A 2.3 CerTest C.*difficile* antigen GDH, Toxin A and Toxin B Positive Control performance e interpretation of the results.

ANNEX 3

KIT 20 CARD *C. DIFFICILE* GDH+TOXIN A+B. Catalogue reference: GX872001VCC

Annex for the following reference:

PRODUCT	REFERENCE
CerTest <i>C. difficile</i> GDH+Toxin A+B Combo Card (20 tests kit with stool vials +PC +NC)	GX872001VCC

Table A 3. Reference for kit 20 combo card of CerTest *C. difficile* GDH+Toxin A+B with positive and negative control.

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
- 20x CerTest <i>C. difficile</i> GDH+Toxin A+B (combo card test format) -20x Stool collection tubes with diluent. -1x CerTest <i>C. difficile</i> antigen GDH Positive Control -1x CerTest <i>C. difficile</i> Toxin A Positive Control -1x CerTest <i>C. difficile</i> Toxin B Positive Control -1x CerTest Negative Control -4x Testing tube -4x Disposable pipette -3x Positive control reagent tubes (red cap) -1x IFU Instructions for use (including Quick guide use).	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Specimen collection container


Table A 3.1 Materials necessities for test performance.

II. Test procedures (samples)


Follow proper infection control practices. Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the **stool collection tube** in order to assure good sample dispersion.
2. Remove **CerTest C. difficile GDH+Toxin A+B** (combo card test) from its sealed bag just before using it.
3. Take the **stool collection tube**, cut the end of the cap (1) and dispense 3 drops in the circular window marked with the letter A (2), 3 drops, using the same tube, in the circular window marked with letter B (3) and 3 drops, using the same tube, in the circular window marked with letter C (4). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.


(2) Add 3 drops in the circular window (A)

(3) Add 3 drops in the circular window (B)





(4) Add 3 drops in the circular window (C)





Wait 10 minutes

III. Interpretation of results

(please, refer to the illustration below)

	GDH	Toxin A	Toxin B		Interpretation of the results
1.	NEGATIVE	NEGATIVE	NEGATIVE		There are no GDH, Toxin A or Toxin B presence. No infection caused by <i>Clostridium difficile</i>. A negative result should be followed up with additional laboratory techniques to confirm the results.
	GREEN	GREEN	GREEN		
2.	POSITIVE	POSITIVE	POSITIVE		There are GDH, Toxin A and Toxin B presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN-RED	GREEN-RED		
3.	POSITIVE	POSITIVE	NEGATIVE		There is GDH and Toxin A presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN-RED	GREEN		
4.	POSITIVE	NEGATIVE	POSITIVE		There is GDH and Toxin B presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN	GREEN-RED		

	GDH	Toxin A	Toxin B		Interpretation of the results
5.	POSITIVE	NEGATIVE	NEGATIVE		There is GDH presence. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN-RED	GREEN	GREEN		
6.	NEGATIVE	POSITIVE	POSITIVE		If this result appears it must be repeat the test using a fresh sample. If results are again positive for Toxin A and B and negative for GDH, the sample should be considered positive for Toxin A and B. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN	GREEN-RED	GREEN-RED		
7.	NEGATIVE	POSITIVE	NEGATIVE		If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin A and negative for GDH, the sample should be considered positive for Toxin A. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	GREEN	GREEN-RED	GREEN		
8.	NEGATIVE	NEGATIVE	POSITIVE		If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin B and negative for GDH, the sample should be considered positive for Toxin B. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line),
	GREEN	GREEN	GREEN-RED		

	GDH	Toxin A	Toxin B		Interpretation of the results
					independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
9.	INVALID				Invalid result either A, B or C, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid.
	Any other result				

Table A 3.2 CerTest C. difficile GDH+Toxin A+B performance interpretation of the results.

Whatever results detailed between 6 and 8, are not very common. A very low percentage of specimens might result negative for GDH but positive for toxins.

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality control

CerTest C. difficile GDH+Toxin A+B (combo card test) presents five controls:

INTERNAL CONTROL: CONTROL LINE

Internal control: internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

EXTERNAL CONTROL: POSITIVE CONTROL SWABS AND NEGATIVE CONTROL SWAB

External positive control: CerTest C. difficile antigen GDH Positive Control, CerTest C. difficile Toxin A Positive Control and CerTest C. difficile Toxin B Positive Control (positive control swab).

External negative control: CerTest Negative Control (negative control swab).

A. External positive controls

It is recommended to use the same device for both external positive controls.

1. CerTest C. difficile antigen GDH Positive Control


CerTest C. difficile antigen GDH Positive Control is an external quality control for **CerTest C. difficile GDH+Toxin A+B** (Strip A). Use **CerTest C. difficile antigen GDH Positive Control** to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C, **CerTest C. difficile antigen GDH Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest C. difficile antigen GDH Positive Control** must remain in the sealed pouch until use. Do not freeze.

1.1. Positive control procedure

Allow tests and CerTest C. difficile antigen GDH Positive Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Take the **positive control reagent tube**, cut the end of the cap.
2. Add 15 drops (1) from the **positive control reagent tube** into the **testing tube** and immediately introduce **CerTest C. difficile antigen GDH Positive Control**.
3. Mix the solution by rotating **CerTest C. difficile antigen GDH Positive Control** forcefully against the side of the **testing tube** at least 15 seconds. Extract as much liquid as possible from the **CerTest C. difficile antigen GDH Positive Control**, squeezing the sides of the **testing tube** as **CerTest C. difficile antigen GDH Positive Control** is withdrawn (2). Discard **CerTest C. difficile antigen GDH Positive Control**.
4. Remove the **CerTest C. difficile GDH+Toxin A+B** (combo card format) from its sealed bag just before using it.
5. Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** in the circular window marked with the letter A (3).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

2. CerTest C. difficile Toxin A Positive Control


CerTest C. difficile Toxin A Positive Control is an external quality control for **CerTest C. difficile GDH+Toxin A+B** (Strip B). Use **CerTest C. difficile Toxin A Positive Control** to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C, **CerTest C. difficile Toxin A Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest C. difficile Toxin A Positive Control** must remain in the sealed pouch until use. Do not freeze.

2.1. Positive control procedure

Allow tests and CerTest C. difficile Toxin A Positive Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Take the **positive control reagent tube**, cut the end of the cap.
2. Add 15 drops (1) from **positive control reagent tube** into the **testing tube** and immediately introduce **CerTest C. difficile Toxin A Positive Control**.
3. Mix the solution by rotating **CerTest C. difficile Toxin A Positive Control** forcefully against the side of the **testing tube** at least 15 seconds. Extract as much liquid as possible from **CerTest C. difficile Toxin A Positive Control**, squeezing the sides of the **testing tube** as **CerTest C. difficile Toxin A Positive Control** is withdrawn (2). Discard **CerTest C. difficile Toxin A Positive Control**.
4. Remove the **CerTest C. difficile GDH+Toxin A+B** (combo card format) from its sealed bag just before using it.
5. Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** in the circular window marked with the letter B (3).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Wait 10 minutes

3. CerTest C. difficile Toxin B Positive Control


CerTest C. difficile Toxin B Positive Control is an external quality control for **CerTest C. difficile GDH+Toxin A+B** (Strip C). Use **CerTest C. difficile Toxin B Positive Control** to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C, **CerTest C. difficile Toxin B Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest C. difficile Toxin B Positive Control** must remain in the sealed pouch until use. Do not freeze.

3.1. Positive control procedure

Allow tests and **CerTest C. difficile Toxin B Positive Control** to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Take the **positive control reagent tube**, cut the end of the cap.
2. Add 15 drops (1) from **positive control reagent tube** into the **testing tube** and immediately introduce the **CerTest C. difficile Toxin B Positive Control**.
3. Mix the solution by rotating the **CerTest C. difficile Toxin B Positive Control** forcefully against the side of the **Testing tube** at least 15 seconds. Extract as much liquid as possible from the **CerTest C. difficile Toxin B Positive Control**, squeezing the sides of the **testing tube** as the **CerTest C. difficile Toxin B Positive Control** is withdrawn (2). Discard **CerTest C. difficile Toxin B Positive Control**.
4. Remove the **CerTest C. difficile GDH+Toxin A+B** (combo card format) from its sealed bag just before using it.
5. Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** in the circular window marked with the letter C (3).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Wait 10 minutes

Positive controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

A.1 Interpretation positive control procedure

	CerTest C. <i>difficile</i> antigen GDH Positive Control	CerTest C. <i>difficile</i> Toxin A Positive Control	CerTest C. <i>difficile</i> Toxin B Positive Control		Interpretation of the results
1.	POSITIVE (One device for three positive controls)				A GREEN line (control line (C)) and a RED line (test line (T)) in the three strips, in the results window, appear at the same time due to presence of <i>C.difficile</i> antigens in CerTest C. <i>difficile</i> antigen GDH Positive Control , in strip A, in the CerTest C. <i>difficile</i> Toxin A Positive Control , in strip B, and in CerTest C. <i>difficile</i> Toxin B Positive Control , in strip C.
	GREEN-RED	GREEN-RED	GREEN-RED		
2.	POSITIVE (One device for CerTest C. <i>difficile</i> antigen GDH Positive Control)				A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window in strip A due to presence of antigen GDH in CerTest C. <i>difficile</i> antigen GDH Positive Control .
	GREEN-RED	-	-		
3.	POSITIVE (One device for CerTest C. <i>difficile</i> Toxin A Positive Control)				A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window in




	CerTest C. <i>difficile</i> antigen GDH Positive Control	CerTest C. <i>difficile</i> Toxin A Positive Control	CerTest C. <i>difficile</i> Toxin B Positive Control		Interpretation of the results
	-	GREEN-RED	-		strip B due to presence of Toxin A in CerTest C.<i>difficile</i> Toxin A Positive Control.
4.	POSITIVE (One device for CerTest C.<i>difficile</i> Toxin B Positive Control)				A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window in strip C due to presence of Toxin B in CerTest C.<i>difficile</i> Toxin B Positive Control.
	-	-	GREEN-RED		
4.	INVALID				Total absence of the control-coloured line (GREEN), regardless the appearance or not of the test line (RED), or only the control (C) coloured line (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new CerTest C. <i>difficile</i> antigen GDH Positive Control, CerTest C. <i>difficile</i> Toxin A Positive Control and CerTest C. <i>difficile</i> Toxin B Positive Control. If the situation still persists, discontinue using the test kit and contact your local distributor. Any result read after more than 10 minutes will be considered invalid.
	Any other result				

Table A 3.3 CerTest C.*difficile* antigen GDH, Toxin A and Toxin B Positive Control performance e interpretation of the results.

B. External negative control. CerTest Negative Control

CerTest Negative Control is an external quality control for **CerTest C. *difficile* GDH+Toxin A+B** (Strip A, B and C). Use **CerTest Negative Control** to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C, **CerTest Negative Control** is stable through the expiration date printed on the sealed pouch. **CerTest Negative Control** must remain in the sealed pouch until use. Do not freeze.

B.1 Negative control procedure


Allow tests and CerTest Negative Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Add 15 drops (1) of diluent from the **stool collection tube** into the **testing tube** and immediately put **CerTest Negative Control**.
2. Mix the solution by rotating **CerTest Negative Control** forcefully against the side of the **testing tube** at least 15 seconds. Extract as much liquid as possible from the **CerTest Negative Control**, squeezing the sides of the **testing tube** as **CerTest Negative Control** is withdrawn (2). Discard the **CerTest Negative Control**.
3. Remove the **CerTest C. difficile GDH+Toxin A+B** (combo card format) from its sealed bag just before using it.
4. Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** into the circular window marked with letter A (3) dispense 3 drops from the **testing tube** into the circular window marked with letter B (4) and dispense 3 drops from the **testing tube** into the circular window marked with letter C (5).
5. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

(1)

Add 15 drops of diluent


(2)

Introduce the **CerTest Negative Control** into the **testing tube**, rotating 15 seconds and extract the liquid.


(2) Add 3 drops in the circular window (A)

(3) Add 3 drops in the circular window (B)


(4) Add 3 drops in the circular window (C)

Control Line Test Line

Control Line Test Line


Control Line Test Line Wait 10 minutes

Negative controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

B.1.1 Interpretation negative control procedure

Negative Control		Interpretation of the results
1.	VALID RESULT	Only A GREEN line (control line (C)) appear across the results window during the test performance.
	GREEN (only control band appear)	

	Negative Control		Interpretation of the results
2.	INVALID	Invalid result either A, B or C, we recommend repeating the assay using the same sample with another test.	Total absence of the control-coloured line (GREEN), regardless the appearance or not of the test line (RED), or present of both lines (RED test and GREEN control lines). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test device and a new CerTest Negative Control. If the situation still persists, discontinue using the test kit and contact your local distributor.
	Any other result		

Table A 3.4 CerTest Negative Control performance interpretation of the results.

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11. Reference/Bibliographie

1. Vaishnavi C. **Clinical spectrum & pathogenesis of Clostridium difficile associated diseases**. Indian J Med Res. 2010;131(4):487–99.
2. Fitzsimmons. **Independent evaluation of a GDH enzyme immunoassay for C. difficile screening**. Biomed Sci. 2011;(October):698–700.
3. Ofori E, Ramai D, Dhawan M, Mustafa F, Gasperino J, Reddy M. **Community-acquired Clostridium difficile: epidemiology, ribotype, risk factors, hospital and intensive care unit outcomes, and current and emerging therapies**. J Hosp Infect [Internet]. 2018;99(4):436–42.
4. Guh AY, Mu Y, Winston LG, Johnston H, Olson D, Farley MM, et al. **Trends in U.S. Burden of Clostridioides difficile Infection and Outcomes**. N Engl J Med. 2021;382(14):1320–30.
5. Wren MWD, Sivapalan M, Kinson R, Shetty NP. **An evaluation of tests for faecal toxin , glutamate dehydrogenase , lactoferrin and toxigenic culture in the diagnostic laboratory**. Br J Biomed Sci. 2009;66(1):1–5.
6. Services U. D of H and H. **Antibiotic resistance threats in the United States**. Centers Dis Control Prev. 2019;1–113.

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