Symbol	EN	FR	ES	IT	DE	PT	HR	SV	NL
	Sufficient for <n> tests	Suffisant pour <n> tests	Suficiente para <n> pruebas	Sufficiente per <n> analisi	Inhalt ausreichend für <n> Prüfungen	Suficiente para <n> testes	Dostatno za <n> testova	Räcker till <n> tester	Voldoende voor <n> tests
	Read instruction for use	Veuillez lire les Instructions d'utilisation	Lea las instrucciones de uso	Leggere le istruzioni per l'uso	Gebrauchsanweisung beachten	Leia as instruções de uso	Pročitati upute za uporabu	Läs bruksanvisningen	Lees de gebruiksaanwijzing
	Use by Date	Date de péremption	Fecha de caducidad	Usare entro	Verwendbar bis	Validade	Upotrijebiti do datuma	Används senast datum	Gebruik voor Datum
	Batch code	Code du lot	Código de lote	Codice lotto	Chargenbezeichnung	Código de lote	Šifra serije	Lotnummer	Lotnummer
	Catalog number	Numéro de catalogue	Número de catálogo	Numero di catalogo	Katalognummer	Número do catálogo	Kataloški broj	Katalognummer	Catalogus nr.
	Caution	Attention	Precaución	Attenzione	Vorsicht	Atenção	Oprez	Försiktighet	Waarschuwing
	Manufacturer	Fabricant	Fabricante	Produttore	Hersteller	Fabricante	Proizvođač	Tillverkare	Fabrikant
	Importer	Importateur	Importador	Importatori	Importeur	Importador	Uvoznik	Importörernas	Importeur
	Authorized representative of the European Community	Représentant autorisé de la Communauté européenne	Representante autorizado en la Comunidad Europea	Rappresentante autorizzato nella Comunità Europea	Bevollmächtigter in der Europäischen Union	Representante autorizado da Comunidade Europeia	Ovlašteni predstavnik Europske zajednice	Behörig företrädare för Europeiska mänskaper	Bevoegde vertegenwoordiger van de Europese Gemeenschap
	In vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Dispositivo médico de diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	In-Vitro-Diagnostikum	Dispositivo médico para diagnóstico in vitro	In vitro dijagnostički medicinski uređaj	Medicinsk utrustning för in vitro-diagnostik	In vitro diagnostisch medisch hulpmiddel
	Temperature limit	Limite de température	Límite de temperatura	Limite di temperatura	Temperaturbegrenzung	Limite de temperatura	Temperaturni režim	Temperatur gräns	Temperatuur grens
	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Nicht zur Wiederverwendung	Não reutilizar	Ne ponovno upotrijebiti	Återanvänd inte	Niet opnieuw gebruiken
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices	Ce produit est conforme aux exigences de la Directive 98/79/EC sur les dispositifs médicaux de diagnostic in vitro	Este producto cumple con los requisitos de la Directiva 98/79/EC sobre dispositivos de diagnóstico in vitro	Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medici diagnostici in vitro	Dieses Produkt entspricht den Anforderungen der europäischen Richtlinie 98/79/EG über In-vitro-Diagnostika	Este produto cumpre os requisitos da Diretiva 98/79/CE sobre dispositivos médicos de diagnóstico in vitro	Ovaj proizvod ispunjava zahtjeve Direktive 98/79/EC o in vitro dijagnostičkim medicinskim uređajima	Denna produkt uppfyller kraven i direktiv 98/79/EG om medicintekniska produkter för in vitro-diagnostik	Dit product voldoet aan de vereisten van de Richtlijn 98/79/EG op in vitro diagnostisch medisch hulpmiddelen

ez-Tracker 

Anti-Etanercept

BIOSYNEX

14 Rue Ambroise Croizat 77183 Croissy Beaubourg France
Tel: +(33)-1-64-62-10-12, Fax: +(33)-1-64-62-09-66
www.biosynex.com

Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon,
Chuncheon-si, Gang-won-do 24398, Republic of Korea
Tel: +(82)-33-243-1400, Fax: +(82)-33-243-9373
www.boditech.co.kr

Obelis s.a.

Bd. Général Wahis 53, 1030 Brussels, Belgium
Tel: +(32)-2-732-59-54, Fax: +(32)-2-732-60-03
mail@obelis.net

ez-Tracker Anti-Etanercept

INTENDED USE

ez-Tracker Anti-Etanercept is a fluorescence immunoassay (FIA) for the quantitative determination of free antibodies to etanercept in human whole blood/serum.

For *in vitro* diagnostic use only.

INTRODUCTION

Etanercept is a dimeric fusion protein composed of the extracellular portion of the human p75 tumor necrosis factor alpha (TNF α) receptor linked to the Fc portion of human immunoglobulin 1 (IgG1) that is used to treat or manage a variety of chronic inflammatory diseases including rheumatoid arthritis (RA), juvenile idiopathic arthritis (JIA) and psoriatic arthritis (PSA), plaque psoriasis and ankylosing spondylitis (AS) by acting as a TNF α inhibitor (TNFi).^[1]

Clinical trials have demonstrated that a low circulating drug level is associated with the lack of clinical response for anti-TNF therapy.^[2,3] The major concern about loss of clinical response is the development of antibodies against anti-TNF drugs. These can lead to bind specifically to drugs and neutralize or eliminate it which decrease its efficacy.^[4]

Therefore, implementation of therapeutic drug monitoring in standard clinical care can have the potential to improve the effectiveness, and safety of treatment with anti-TNF drugs.^[5,6]

ez-Tracker Anti-Etanercept is used for a quantitative detection of free antibodies to etanercept in serum and whole blood. This test could provide a relevant information about appropriate therapeutic treatment.

PRINCIPLE

ez-Tracker Anti-Etanercept uses a bridging immunoassay for detection of antibodies to etanercept. Anti-**ez-Tracker** antibodies from sample are bridged with fluorescence and biotin labeled etanercept, forming immune complexes and migrate onto nitrocellulose matrix to be captured by the other immobilized streptavidin on a test strip.

More antibodies to etanercept in the sample will form more immune complexes which lead to stronger fluorescence signal, which is processed by the instrument for **ez-Tracker** tests to show free anti-**ez-Tracker** concentration in the sample.

COMPONENT

ez-Tracker Anti-Etanercept consists of 'cartridges'.

• Each sealed aluminum pouch contains two cartridges.

• Each cartridge packaged in an aluminum pouch has three components including a cartridge part, a detector part, and a diluent part.

• The cartridge part contains the membrane called a test strip which has streptavidin at the test line, and chicken IgY at the control line.

• The detector part has 2 granules containing etanercept-fluorescence conjugate, anti-chicken IgY-fluorescence conjugate, etanercept-biotin conjugate, bovine serum albumin (BSA) as a stabilizer in tris-hydrochloride buffer.

• The diluent part contains tween 20 as a detergent and sodium azide as a preservative in phosphate buffered saline (PBS).

WARNING AND PRECAUTIONS

• For *in vitro* diagnostic use only.

• Follow the instructions and procedures described in this 'Instructions for use'.

• Use only fresh samples and avoid direct sunlight.

• Lot numbers of all the test components (cartridge and ID chip) must match each other.

• Do not interchange the test components between the different lots or use the test components after the expiration date, either of which might yield incorrect test result(s).

• Do not reuse cartridges. A cartridge should be used for testing one sample only.

• The cartridge should remain sealed in its original pouch until just before use. Do not use a cartridge, if the pouch is damaged or has already been opened.

• Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used.

• If test components and/or sample are stored in refrigerator, then allow cartridge and sample to be at room temperature for approximately 30 minutes before use.

• The instrument for **ez-Tracker** tests may generate slight vibration during use.

• Used cartridges, C-tips and pipette tips should be handled carefully and discarded by an appropriate method in accordance with relevant local regulations.

• The cartridge contains sodium azide (NaN₃), and they may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure. Avoid contact with skin, eyes, and clothing. In case of contact, rinse immediately with running water.

• No biohazard interference was observed in **ez-Tracker Anti-Etanercept** when biotin concentration in the sample was below 500 ng/mL. If a patient has been taking biotin at dosage of more than 0.03 mg a day, it

is recommended to test again 24 hours after discontinuation of biotin intake.

• **ez-Tracker Anti-Etanercept** will provide accurate and reliable results subject to the below conditions.

- **ez-Tracker Anti-Etanercept** should be used only in conjunction with the instrument for **ez-Tracker** tests.

- Have to use recommended anticoagulant.

Recommended anticoagulant	
K ₂ EDTA, Na ₂ EDTA, Sodium Citrate, Lithium heparin, Sodium heparin	

• **C-tip should be used when the following conditions are met.**

- C-tip provided with the kit is recommended to obtain correct test result.

- Whole blood should be immediately tested after collection.

- Do not perform a test with C-tip on General Mode. It might cause an erroneous result.

- Excess whole blood around the C-tip should be wiped off.

- In order to avoid cross-contamination, please do not re-use C-tip for multiple samples.

- Cartridge should be inserted and positioned in the cartridge holder prior to the blood sample collection.

- While collecting blood, be careful not to create air bubbles in the C-tip.

LIMITATIONS OF THE TEST SYSTEM

• The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.

• The test may yield false negative result(s) due to the non-responsiveness of the anti-drug antibody which is the most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the drug with time and/or temperature may also cause false negative result as it makes the drug unrecognizable by the antibodies.

• Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.

• Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.

STORAGE AND STABILITY

Component	Storage condition		Note
	Storage Temperature	Shelf life	
Cartridge	2 - 30 °C	20 months 1 months	Unopened Resealed

• Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal.

MATERIALS SUPPLIED

REF ETE 003-24

Components of **ez-Tracker Anti-Etanercept**

Cartridge	CART	24
Pipette tip (Zipper bag)	P-T	24
C-tip (Zipper bag)	C-T	24
ID chip	ID-C	1
Spare cartridge zipper bag	SCZB	1
Instructions for use	IFU	1

MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately with **ez-Tracker Anti-Etanercept**.

Please contact our sales division for more information.

• **ez-Track[®]**

• **ez-Tracker Anti-Etanercept Control**

• **ez-Tracker Anti-Etanercept Calibrator**

REF ETE 001

REF ETE 003-C

REF ETE 003-CAL

SAMPLE COLLECTION AND PROCESSING

The sample type for **ez-Tracker Anti-Etanercept** is human whole blood/serum.

• It is recommended to test the sample within 24 hours after collection.

• The serum should be separated from the clot by centrifugation within 3 hours after the collection of whole blood.

• The serum may be stored for a week at 2 - 8 °C prior to being tested. If testing will be delayed more than a week, samples should be frozen at -20 °C.

• The serum stored frozen at -20 °C for 1 month showed no performance difference.

• However, the whole blood sample should not be kept in a freezer in any case.

• As a repeated freeze-thaw cycle may affect the test result, do not refreeze previously frozen samples.

• Collection of whole blood sample using C-tip

① Hold the C-tip horizontally and touch the surface of the blood drop with the tip of the C-tip.

② Capillary action will automatically draw the blood sample to C-tip and stop.

③ Wipe off any excess blood around the tip.

④ Double-check if whole blood is filled accurately in the C-tip and instrument for **ez-Tracker** tests is ready for a test on the 'C-tip' mode.

TEST SETUP

• Check the components of the **ez-Tracker Anti-Etanercept** as described below: Cartridges, pipette tips, C-tips, an ID chip, a spare cartridge zipper bag and an instructions for use.

• Ensure that the lot number of the cartridge matches that of the ID chip.

- the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for ez-Tracker tests.
- Empty the tip box.
- Insert the ID chip into the 'ID chip port'.

※ **Please refer to the instrument for ez-Tracker tests operation manual for complete information and operating instructions.**

TEST PROCEDURE

General mode

- 1) Insert a cartridge into the cartridge holder.
- 2) Insert a tip into the tip hole of the cartridge.
- 3) Select the 'General mode' in the instrument for ez-Tracker tests.
- 4) Take 100 μ L of the sample (whole blood/serum/control) using a pipette and dispense it into the sample well of the cartridge.
- 5) Tap the 'Start' button on the screen.
- 6) The test result will be displayed on the screen after 10 minutes.

C-tip mode

- 1) Insert a cartridge into the cartridge holder.
- 2) Take 30 μ L of whole blood using a C-tip.
- 3) Insert the whole blood-filled C-tip into the tip hole of the cartridge.
- 4) Select the 'C-tip mode' in the instrument for ez-Tracker tests.
- 5) Tap the 'Start' button on the screen.
- 6) The test result will be displayed on the screen after 10 minutes.

INTERPRETATION OF TEST RESULT

- The instrument for ez-Tracker tests calculates the test result automatically and displays Anti-Etanercept concentration of the test sample in terms ng/mL.
- Working range: 10 - 100 ng/mL

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand with ez-Tracker Anti-Etanercept. For more information regarding obtaining the control materials, contact Biosynex's Sales Division for assistance.
(Please refer to the instructions for use of control material.)

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

- Limit of Blank (LOB)	3.06 ng/mL
- Limit of Detection (LOD)	4.22 ng/mL
- Limit of Quantitation (LOQ)	10.00 ng/mL

High-dose Hook Effect

There is no high-dose hook effect at anti-etanercept concentration up to 30,000 ng/mL.

Analytical specificity

Cross-reactivity

ez-Tracker Anti-Etanercept test results did not show any significant cross-reactivity with these biomolecules.

Cross-reactants	Concentration
Infliximab	100 μ g/mL
Adalimumab	100 μ g/mL
Vedolizumab	100 μ g/mL
Golimumab	100 μ g/mL
Anti-infliximab	1,000 ng/mL
Anti-adalimumab	1,000 ng/mL
Anti-vedolizumab	1,000 ng/mL
Anti-golimumab	1,000 ng/mL

Interference

Interferents listed in the following table were added to the test sample at the concentration mentioned below. ez-Tracker Anti-Etanercept test results did not show any significant interference with these materials.

Interferents	Concentration
Hemoglobin	1,000 mg/dL
Bilirubin	40 mg/dL
Triglyceride	1,500 mg/dL
Rheumatoid factor	80 IU/mL
Human serum albumin	12 g/dL

Precision

Single-site study

Repeatability (within-run precision)
within-laboratory precision (total precision)
Lot to lot precision

3 Lots of ez-Tracker Anti-Etanercept were tested for 20 days. Each standard material was tested 2 times per day. For each test, each material was duplicated.

Multi-site study

Reproducibility

1 Lot of ez-Tracker Anti-Etanercept was tested for 5 days in 3 different sites (1 person per site, 1 instrument per site). Each standard material was tested 1 time per and 5 replicates per day.

Anti-Etanercept	Repeatability		Total precision	
[ng/mL]	AVG	CV (%)	AVG	CV (%)
15	15.19	6.6	15.09	6.8
30	29.08	8.1	29.23	7.6
60	58.52	6.9	58.99	6.4
Anti-Etanercept	Lot to lot precision		Reproducibility	
[ng/mL]	AVG	CV (%)	AVG	CV (%)
15	14.98	6.7	15.03	6.2
30	29.14	8.0	29.51	7.8
60	58.83	6.3	60.69	5.4


Accuracy

The accuracy was confirmed by testing with 3 different lots of ez-Tracker Anti-Etanercept. The tests were repeated 10 times at each concentration of the control standard.

Anti-Etanercept [ng/mL]	Lot 1	Lot 2	Lot 3	AVG	Recovery (%)
22.14	23.13	21.90	22.31	22.44	101.4
35.83	36.39	36.72	35.99	36.37	101.5
57.50	57.56	56.97	57.40	57.31	99.7
71.43	73.13	75.75	69.54	72.80	101.9
86.67	91.19	87.60	90.02	89.60	103.4

Comparability

Anti-Etanercept concentration of 74 clinical samples were quantified independently with ez-Tracker Anti-Etanercept and comparator A as per prescribed test procedures. Test results were compared, and their comparability was investigated with linear regression and correlation coefficient (R). The regression equation and correlation coefficient are as follows.

REFERENCES

1. Meenu Wadhwa, Chris Bird, Paula Dilger, Peter Rigby et al. Establishment of the first WHO International Standard for etanercept, a TNF receptor II C fusion protein: Report of an international collaborative study. *J Immunol Methods*. 2017 Aug;447:14-22.
2. Bartelds GM, Kriekamp CT, Nurmohamed MT, et al. Development of anti-tumor necrosis factor against adalimumab and association with disease activity and treatment failure during long-term follow-up. *JAMA* 2011;305:1460-8.
3. Wolinski GJ, Vis M, Lems W, et al. Development of anti-infliximab antibodies and relationship to clinical response in patients with rheumatoid arthritis. *Arthritis Rheum* 2006;54:711-15.
4. A. Nassar-Sheikh Rashid, D. Schonenberg-Meinema, S. C. Bergkamp et al. Therapeutic drug monitoring of anti-TNF drugs: an overview of applicability in daily clinical practice in the era of treatment with biologics in juvenile idiopathic arthritis (JIA) *Pediatric Rheumatology* volume 19, Article number: 59 (2021)
5. Martin-Lopez M, Carmona L, Balsa A, Calvo-Alen J, Santamari R, Tornero J, et al. Serum drug levels of biologic agents in the management of rheumatoid arthritis and spondyloarthritis: a systematic review. *Rheumatol Int* 2018;38:975-83.
6. Kriekamp CT, Nair SC, Nurmohamed MT, van Dongen CJ, Lems WF, Lafere FP, et al. Personalised treatment using serum drug levels of adalimumab in patients with rheumatoid arthritis: an evaluation of costs and effects. *Ann Rheum Dis* 2015;74:361-8.
7. Mullemann D, Meric JC, Paintaud G, Ducourau E, Magdelaine Beuzelin C, Valat JP, et al. Infliximab concentration monitoring improves the control of disease activity in rheumatoid arthritis. *Arthritis Res Ther* 2009;11:R178.
8. IG Sanchez-Hernandez, N Rebollo, F Munoz et al. Therapeutic drug monitoring of tumour necrosis factor inhibitors in the management of chronic inflammatory diseases. *Annals of Clinical Biochemistry* 2019, Vol. 56(1) 28-41.
9. Papamichael K, Cheifetz AS. Use of anti-TNF drug levels to optimise patient management. *Frontline Gastroenterology* 2016;0:1-12.

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
RÉFÉRENCES

1. Meenu Wadhwa, Chris Bird, Paula Dilger, Peter Rigby et al. Establishment of the first WHO International Standard for etanercept, a TNF receptor II Fc fusion protein: Report of an international collaborative study. *J Immunol Methods*. 2017; 447:14-22.
2. Bartelds GW, Kriekkaert CJ, Nurmohamed MT, et al. Development of anti-drug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. *JAMA* 2011;305:1460-8.
3. Wolbink GJ, Vis M, Lems W, et al. Development of anti-infliximab antibodies and relationship to clinical response in patients with rheumatoid arthritis. *Arthritis Rheum* 2006;54:711-15.
4. A. Nassar-Sheikh Rashid, D. Schonenberg-Meinema, S. C. Bengkamp et al. Therapeutic drug monitoring of anti-TNF drugs: an overview of applicability in daily clinical practice in the era of treatment with biologics in juvenile idiopathic arthritis (JIA) *Pediatric Rheumatology* volume 19, Article number: 59 (2021)
5. Martin-Lopez M, Carmona L, Balsa A, Calvo-Alen J, Sanmarti R, Tornero J, et al. Serum drug levels of biologic agents in the management of rheumatoid arthritis and spondyloarthritis: a systematic review. *Rheumatol Int* 2018;38:975-83
6. Kriekkaert CJ, Nair SC, Nurmohamed MT, van Dongen CJ, Lems WF, Lefebvre FP, et al. Personalised treatment using serum drug levels of adalimumab in patients with rheumatoid arthritis: an evaluation of costs and effects. *Ann Rheum Dis* 2015;74:361-8.
7. Mullemann D, Meric JC, Paintaud G, Ducourau A, Magdeleine Beuzelin C, Valat JP, et al. Infliximab concentration monitoring improves the control of disease activity in rheumatoid arthritis. *Arthritis Res Ther* 2009;11:R78.
8. IG Sanchez-Hernandez, N. Rebollo, F Munoz et al. Therapeutic drug monitoring of tumour necrosis factor inhibitors in the management of chronic inflammatory diseases. *Annals of Clinical Biochemistry* 2019; Vol. 56(1) 28-41.
9. Papamichael K, Cheifetz AS. Use of anti-TNF drug levels to optimise patient management. *Frontline Gastroenterology* 2016;0:1-12.

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
REFERENCIAS

1. Meenu Wadhwa, Chris Bird, Paula Dilger, Peter Rigby et al. Establishment of the first WHO International Standard for etanercept, a TNF receptor II Fc fusion protein: Report of an international collaborative study. *J Immunol Methods*. 2017 Aug;447:14-22.
2. Bartelds GW, Kriekraert CL, Nurmohamed MT, et al. Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. *JAMA*. 2011;305:1469-8.
3. Wolinski GJ, Vis M, Lems W, et al. Development of anti infliximab antibodies and relationship to clinical response in patients with rheumatoid arthritis. *Arthritis Rheum* 2006;54:711-15.
4. A. Nassar-Sheikh Rashid, D. Schonenberg-Meinema, S. C. Bergkamp et al. Therapeutic drug monitoring of anti-TNF drugs: an overview of applicability in daily clinical practice in the era of treatment with biologics in juvenile idiopathic arthritis (JIA). *Pediatric Rheumatology volume 19, Article number: 59* (2021)
5. Martin-Lopez M, Carmona L, Balsa A, Calvo-Alen J, Sanmartín R, Tornero J, et al. Serum drug levels of biologic agents in the management of rheumatoid arthritis and spondyloarthritis: a systematic review. *Rheumatol Int* 2018;38:975-83.
6. Kriekraert CL, Nair SC, Nurmohamed MT, van Dongen CJ, Lems WF, Lafeyre FP, et al. Personalised treatment using serum drug levels of adalimumab in patients with rheumatoid arthritis: an evaluation of costs and effects. *Ann Rheum Dis* 2015;74:361-8.
7. Mullemann D, Meric JC, Paintaud G, Ducourau E, Magdelaine Beuzelin C, Valat JP, et al. Infliximab concentration monitoring improves the control of disease activity in rheumatoid arthritis. *Arthritis Res Ther* 2009;11:R178.
8. IG Sanchez-Hernandez, N Rebollo, F Munoz et al. Therapeutic drug monitoring of tumour necrosis factor inhibitors in the management of chronic inflammatory diseases. *Annals of Clinical Biochemistry* 2019; Vol. 56(1) 28-41.
9. Papamichael K, Cheifetz AS. Use of anti-TNF drug levels to optimise patient management. *Frontline Gastroenterology* 2016;0:1-12.

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
BIBLIOGRAFIA

- Meenu Wadhwa, Chris Bird, Paula Dilger, Peter Rigby et al. Establishment of the first WHO International Standard for etanercept, a TNF receptor 1 Fc fusion protein: Report of an international collaborative study. *J Immunol Methods*. 2017 Aug;447:14-22.
- Bartelds GM, Kriekciak CT, Nurmohamed MT, et al. Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. *JAMA*. 2011;305:1460-8.
- Wolbink GJ, Vis M, Lems W, et al. Development of anti-infliximab antibodies and relationship to clinical response in patients with rheumatoid arthritis. *Arthritis Rheum* 2006;54:711-15.
- A. Nassar-Sheikh Rashid, D. Schonenberg-Meinema, S. C. Berkamp et al. Therapeutic drug monitoring of anti-TNF drugs: an overview of applicability in daily clinical practice in the era of treatment with biologics in juvenile idiopathic arthritis (JIA). *Pediatric Rheumatology* volume 19, Article number: 59 (2012)
- Martin-Lopez M, Carmona L, Balsa A, Calvo-Alen J, Sanmartín R, Tornero J, et al. Serum drug levels of biologic agents in the management of rheumatoid arthritis and spondyloarthritis: a systematic review. *Rheumatol Int* 2018;38:975-83.
- Kriekciak CT, Nair SC, Nurmohamed MT, van Dongen CJ, Lems WF, Lefebvre FP, et al. Personalised treatment using serum drug levels of adalimumab in patients with rheumatoid arthritis: an evaluation of costs and effects. *Ann Rheum Dis* 2015;74:361-8.
- Mulleman D, Meric JC, Paintaud G, Ducourau A, Magdelaine Beuzelin C, Valat JP, et al. Infliximab concentration monitoring improves the control of disease activity in rheumatoid arthritis. *Arthritis Res Ther* 2009;11:R178.
- JG Sanchez-Hernandez, N Rebollo, F Munoz et al. Therapeutic drug monitoring of tumour necrosis factor inhibitors in the management of chronic inflammatory diseases. *Annals of Clinical Biochemistry* 2019, Vol. 56(1) 28-41.
- Papamichael K, Cheifetz AS. Use of anti-TNF drug levels to optimise patient management. *Frontline Gastroenterology* 2016;1-12.

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Kreuzreaktive Stoffe	Konzentration
Infliximab	100 µg/ml
Adalimumab	100 µg/ml
Vedolizumab	100 µg/ml
Golimumab	100 µg/ml
Anti-Infliximab	1,000 ng/ml
Anti-Adalimumab	1,000 ng/ml
Anti-Vedolizumab	1,000 ng/ml
Anti-Golimumab	1,000 ng/ml

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
QUELLENANGABEN

- Meenu Wadhwa, Chris Bird, Paula Dilger, Peter Rigby et al. Establishment of the first WHO International Standard for etanercept, a TNF receptor II Fc fusion protein: Report of an international collaborative study. *J Immunol Methods*. 2017; 447:14-22.
- Bartelds GM, Kriekciak CT, Nurmohamed MT, et al. Development of anti-tumor necrosis factor (TNF) antibodies and association with disease activity and treatment failure during long-term follow-up. *JAMA*. 2011; 305:1460-8.
- Wolbink GJ, Vis M, Lems W, et al. Development of anti-TNF antibodies and relationship to clinical response in patients with rheumatoid arthritis. *Arthritis Rheum*. 2006; 54:711-15.
- A. Nassar-Sheikh Rashid, D. Schonenberg-Meimann, S. C. Bergkamp et al. Therapeutic drug monitoring of anti-TNF drugs: an overview of applicability in daily clinical practice in the era of treatment with biologics in juvenile idiopathic arthritis (JIA) *Pediatric Rheumatology* volume 19, Article number: 59 (2021)
- Martín-López M, Carmona I, Balsa A, Calvo-Alén J, Sanmartín R, Tornero J, et al. Serum drug levels of biologic agents in the management of rheumatoid arthritis and spondyloarthritis: a systematic review. *Rheumatol Int*. 2018; 38:975-83
- Krickack CT, Nair SC, Nurmohamed MT, van Dongen CJ, Lems WF, Laferte PP, et al. Personalised treatment using serum drug levels of adalimumab in patients with rheumatoid arthritis: an evaluation of costs and effects. *Ann Rheum Dis*. 2015; 74:361-8.
- Mulleman D, Meric JC, Pintauro G, Ducourau E, Magdelaine Beuzelin C, Valat JP, et al. Infliximab concentration monitoring improves the control of disease activity in rheumatoid arthritis. *Arthritis Res Ther*. 2009; 11:R178.
- JG Sanchez-Hernandez, N Rebollo, F Munoz et al. Therapeutic drug monitoring of tumour necrosis factor inhibitors in the management of chronic inflammatory diseases. *Annals of Clinical Biochemistry* 2019, Vol. 56(1) 28-41
- Papamichal K, Cheifetz AS. Use of anti-TNF drug levels to optimise patient management. *Frontline Gastroenterology* 2016; 0:1-12.

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
REFERÊNCIAS

1. Meenu Wadhwa, Chris Bird, Paula Dilger, Peter Rigby et al. Establishment of the first WHO International Standard for etanercept, a TNF receptor II Fc fusion protein: Report of an international collaborative study. J Immunol Methods. 2017; 447:14-22.
2. Bartelds GM, Kriekamp AT, Nurmohamed MT, et al. Development of anti-tumor antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. JAMA. 2011; 305:1460-8.
3. Wolinski GJ, Vis M, Lems W, et al. Development of anti-inflammasome antibodies and relationship to clinical response in patients with rheumatoid arthritis. Arthritis Rheum 2006; 54:711-15.
4. A. Nassar-Sheikh Rashid, D. Schonenberg-Meinema, S. C. Bergkamp et al. Therapeutic drug monitoring of anti-TNF drugs: an overview of applicability in daily clinical practice in the era of treatment with biologics in juvenile idiopathic arthritis (JIA) Pediatric Rheumatology volume 19, Article number: 59 (2021)
5. Martin-Lopez M, Carmona L, Balsa A, Calvo-Alen J, Sanmartín R, Tornero J, et al. Serum drug levels of biologic agents in the management of rheumatoid arthritis and spondyloarthritis: a systematic review. Rheumatol Int 2018; 38:975-83.
6. Kriekamp AT, Nair SC, Nurmohamed MT, van Dongen CJ, Lems WF, Lafeyre FP, et al. Personalised treatment using serum drug levels of adalimumab in patients with rheumatoid arthritis: an evaluation of costs and effects. Ann Rheum Dis 2015; 74:361-8.
7. Mullemann D, Meric JG, Pintauro D, Ducourau E, Magdelaine Beuzelin C, Valat JP, et al. Infliximab concentration monitoring improves the control of disease activity in rheumatoid arthritis. Arthritis Res Ther 2009; 11:R178.
8. IG Sanchez-Hernandez, N Rebollo, F Munoz et al. Therapeutic drug monitoring of tumour necrosis factor inhibitors in the management of chronic inflammatory diseases. Annals of Clinical Biochemistry 2019, Vol. 56(1) 28-41.
9. Papamichael K, Cheifetz AS. Use of anti-TNF drug levels to optimise patient management. Frontline Gastroenterology 2016; 0:1-12.

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
REFERENCE

1. Meenu Wadhwa, Chris Bird, Paula Dilger, Peter Rigby et al. Establishment of the first WHO International Standard for etanercept, a TNF receptor II Fc fusion protein: Report of an international collaborative study. J Immunol Methods. 2017 Aug;447:14-22.
2. Bartelds GM, Kriekciak CT, Nurmohamed MT, et al. Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. JAMA. 2011;305:1460-8.
3. Wolbink GJ, Vis M, Lems W, et al. Development of anti-infliximab antibodies and relationship to clinical response in patients with rheumatoid arthritis. Arthritis Rheum 2006;54:711-15.
4. A. Nassar-Sheikh Rashid, D. Schonenberg-Meinema, S. C. Beringkamp et al. Therapeutic drug monitoring of anti-TNF drugs: an overview of applicability in daily clinical practice in the era of treatment with biologics in juvenile idiopathic arthritis (JIA) Pediatric Rheumatology volume 19, Article number: 59 (2021)
5. Martin-Lopez M, Carmona L, Balsa A, Calvo-Alen J, Sanmartín R, Tornero J, et al. Serum drug levels of biologic agents in the management of rheumatoid arthritis and spondyloarthritis: a systematic review. Rheumatol Int 2018;38:975-83.
6. Kriekciak CT, Nair SC, Nurmohamed MT, van Dongen CJ, Lems WF, Lafeyre FP, et al. Personalised treatment using serum drug levels of adalimumab in patients with rheumatoid arthritis: an evaluation of costs and effects. Ann Rheum Dis 2015;74:361-8.
7. Mulleman D, Meric JC, Paintaud G, Ducourau E, Magdeleine Beuzelin C, Valat JP, et al. Infliximab concentration monitoring improves the control of disease activity in rheumatoid arthritis. Arthritis Res Ther 2009;11:R178.
8. IG Sanchez-Hernandez, N. Rebollo, F Munoz et al. Therapeutic drug monitoring of tumour necrosis factor inhibitors in the management of chronic inflammatory diseases. Annals of Clinical Biochemistry 2019, Vol. 56(1) 28-41
9. Papamichael K, Cheifetz AS. Use of anti-TNF drug levels to optimise patient management. Frontline Gastroenterology 2016;10:1-12.

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
REFERENSER

1. Meenu Wadhwa, Chris Bird, Paula Dilger, Peter Rigby et al. Establishment of the first WHO International Standard for etanercept, a TNF receptor II Fc fusion protein: Report of an international collaborative study. J Immunol Methods. 2017 Aug;447:14-22.
2. Bartelds GM, Kriekciak CT, Nurmohamed MT, et al. Development of anti-TNF antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. JAMA. 2011;305:1460-8.
3. Wolbink GJ, Vis M, Lems W, et al. Development of anti-TNF antibodies and relationship to clinical response in patients with rheumatoid arthritis. Arthritis Rheum 2006;54:711-15.
4. A. Nassar-Sheikh Rashid, D. Schonenberg-Meinema, S. C. Bergkamp et al. Therapeutic drug monitoring of anti-TNF drugs: an overview of applicability in daily clinical practice in the era of treatment with biologics in juvenile idiopathic arthritis (JIA) Pediatric Rheumatology volume 19, Article number: 59 (2021)
5. Martin-Lopez M, Carmona L, Balsa A, Calvo-Alen J, Sanmartín R, Tornero J, et al. Serum drug levels of biologic agents in the management of rheumatoid arthritis and spondyloarthritis: a systematic review. Rheumatol Int 2018;38:975-983
6. Kriekciak CT, Nair SC, Nurmohamed MT, van Dongen CJ, Lems WF, Lafeyre FP, et al. Personalised treatment using serum drug levels of adalimumab in patients with rheumatoid arthritis: an evaluation of costs and effects. Ann Rheum Dis 2015;74:361-8.
7. Mullemann D, Meric JC, Paintaud G, Ducourau E, Magdeleine Beuzelin C, Valat JP, et al. Infliximab concentration monitoring improves the control of disease activity in rheumatoid arthritis. Arthritis Res Ther 2009;11:R178.
8. IG Sanchez-Hernandez, N Rebollo, F Munoz et al. Therapeutic drug monitoring of tumour necrosis factor inhibitors in the management of chronic inflammatory diseases. Annals of Clinical Biochemistry 2010, Vol. 56(1) 28-41
9. Papamichael K, Cheifetz AS. Use of anti-TNF drug levels to optimise patient management. Frontline Gastroenterology 2016;10:1-12.

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Geteste biologische	Concentratie
Infliximab	100 µg/mL
Adalimumab	100 µg/mL
Vedolizumab	100 µg/mL
Golimumab	100 µg/mL
Anti-Infliximab	1,000 ng/mL
Anti-Adalimumab	1,000 ng/mL
Anti-Vedolizumab	1,000 ng/mL
Anti-Golimumab	1,000 ng/mL

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
REFERENTIES

1. Meenu Wadhwa, Chris Bird, Paula Dilger, Peter Rigby et al. Establishment of the first WHO International Standard for etanercept, a TNF receptor II Fc fusion protein: Report of an international collaborative study. J Immunol Methods. 2017; 447:14-22.
2. Bartelds GW, Kriekamp AT, Nurmohamed MT, et al. Development of anti-TNF antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. JAMA. 2011; 305:1460-8.
3. Wolinski GJ, Vls M, Lems W, et al. Development of anti-infliximab antibodies and relationship to clinical response in patients with rheumatoid arthritis. Arthritis Rheum 2006; 54:711-15.
4. A. Nasar-Sheikh Rashid, D. Schonenberg-Meinema, S. C. Bergkamp et al. Therapeutic drug monitoring of anti-TNF drugs: an overview of applicability in daily clinical practice in the era of treatment with biologics in juvenile idiopathic arthritis (JIA) Pediatric Rheumatology volume 19, Article number: 59 (2021)
5. Martin-Lopez M, Carmona L, Balsa A, Calvo-Alen J, Samartí R, Tornero J, et al. Serum drug levels of biologic agents in the management of rheumatoid arthritis and spondyloarthritis: a systematic review. Rheumatol Int 2018; 38:975-83.
6. Kriekamp AT, Nair SC, Nurmohamed MT, van Dongen CJ, Lems WF, Lafere FP, et al. Personalised treatment using serum drug levels of adalimumab in patients with rheumatoid arthritis: an evaluation of costs and effects. Ann Rheum Dis 2015; 74:361-8.
7. Mullemann D, Meric JC, Pintauro G, Ducourau E, Magdelaine Beuzelin C, Valat JP, et al. Infliximab concentration monitoring improves the control of disease activity in rheumatoid arthritis. Arthritis Res Ther 2009; 11:R178.
8. IG Sanchez-Hernandez, N Rebollo, F Munoz et al. Therapeutic drug monitoring of tumour necrosis factor inhibitors in the management of chronic inflammatory diseases. Annals of Clinical Biochemistry 2019, Vol. 56(1) 28-41.
9. Papamichael K, Cheifetz AS. Use of anti-TNF drug levels to optimise patient management. Frontline Gastroenterology 2016; 0:1-12.