



For assay validation


A complete range of ready to use, internal auto-immune Quality Control sera, multiparametric, compatible with several techniques (multiplex FIDIS™, ELISA, DOT, ...) and CE marked




CONNECTIVITIS




AUTOIMMUNE
THYROIDITIS



RHEUMATOID
ARTHRITIS




ANTI-
PHOSPHOLIPIDS
SYNDROME




VASCULITIS



INFLAMMATORY
BOWEL
DISEASE (IBD)




CELIAC
DISEASE



PERNICIOUS
ANEMIA



DIABETES
MELLITUS TYPE 1



AUTOIMMUNE
LIVER DISEASES

IMMUNO-TROL sera are multiparametric Internal Quality Controls (IQC) for the routine validation of Autoimmune antibody assay techniques including FIDIS™ (xMAP® technology), ELISA and DOT.

Assigned ranges

- Liquid control, ready to use
- Standardised controls
- Identical use and reactivity as patients samples

Convenient

- Universal multi-parametric controls
- Stable material suitable for intra and inter assay validation

Advantages

- Monitor potential drift of analytical systems
- Quality assurance for assay validation



		dsDNA	SS-A 60 Kda	TRIM 21 (SS-A 52 Kda)	SS-B	Sm	Sm/RNP	U1-RNP	Scl-70	Centromere protein B	Jo-1	Ribosomal proteins	Cardiolipin IgG	β2 GPI IgG	Prothrombin IgG	Phospholipids IgG	Cardiolipin IgM	β2 GPI IgM	Prothrombin IgM	Phospholipids IgM	Animal and human RF IgM	CCP IgG	TPO	TG	R-TSH	Mitochondria	MPO	PR3	GBM	PCA	IF	Transglutaminase IgG	Gliadin IgG	Deamidated Gliadin IgG	Transglutaminase IgA	Gliadin IgA	Deamidated Gliadin IgA	ASCA IgA	ASCA IgG	GAD	IA-2	ZnT8		
CONNECTIVE 1	HM060																																											
CONNECTIVE 2	HM061																																											
APS IgG	HM062																																											
APS IgM	HM063																																											
RHEUMA	HM064																																											
THYRO	HM065																																											
CYTO	HM066																																											
VASCULITIS	HM067																																											
GASTRITIS	HM068																																											
CELIAC G	HM069																																											
CELIAC A	HM070																																											
ASCA	HM071																																											
DIABETE	HM072																																											

Reference

Designation

Parameters

Packaging

IMMUNO-TROL

Multiparametric Internal Autoimmune Quality Controls

HM 059-2X1	IMMUNO-TROL NEG	Auto-immune negative Internal Quality Control	2 x 1 mL
HM 060-2X05	IMMUNO-TROL CONNECTIVE 1	dsDNA, SS-A 60kDa, TRIM 21, SS-B, Sm, Sm/RNP & U1-RNP	2 x 0,5 mL
HM 061-2X05	IMMUNO-TROL CONNECTIVE 2	dsDNA, Scl-70, Centromere protein B, Jo-1 & Ribosomal proteins	2 x 0,5 mL
HM 062-2X05	IMMUNO-TROL APS IgG	Cardiolipin IgG, β2 GPI IgG, Prothrombin IgG & Phospholipids IgG	2 x 0,5 mL
HM 063-2X05	IMMUNO-TROL APS IgM	Cardiolipin IgM, β2 GPI IgM, Prothrombin IgM & Phospholipids IgM	2 x 0,5 mL
HM 064-2X05	IMMUNO-TROL RHEUMA	Animal and Human RF IgM & CCP IgG	2 x 0,5 mL
HM 065-2X05	IMMUNO-TROL THYRO	TPO, TG & R-TSH	2 x 0,5 mL
HM 066-2X05	IMMUNO-TROL CYTO	Ribosomal proteins, Mitochondria	2 x 0,5 mL
HM 067-2X05	IMMUNO-TROL VASCULITIS	MPO, PR3 & GBM	2 x 0,5 mL
HM 068-2X05	IMMUNO-TROL GASTRITIS	PCA & IF	2 x 0,5 mL
HM 069-2X05	IMMUNO-TROL CELIAC G	Transglutaminase IgG, Gliadin IgG & Deamidated Gliadin IgG	2 x 0,5 mL
HM 070-2X05	IMMUNO-TROL CELIAC A	Transglutaminase IgA, Gliadin IgA & Deamidated Gliadin IgA	2 x 0,5 mL
HM 071-2X05	IMMUNO-TROL ASCA	Anti-Saccharomyces Cerevisiae Antibodies IgA & IgG	2 x 0,5 mL
HM 072-2X1	IMMUNO-TROL DIABETE	GAD, IA-2 & ZnT8	2 x 1 mL

Theradiag - 14 rue Ambroise Croizat, CS 90136 Croissy-Beaubourg

F-77435 MARNE LA VALLEE CEDEX 2

Tél. 00 33 (0) 1 64 62 10 12 - Fax 00 33 (0) 1 64 62 09 66

info@theradiag.com - www.theradiag.com

Read carefully the instruction for use of the product insert