

ez-Track¹

***Rapid Fluorescent Assay
for the TDM* of Biologics***

Simple

- 1 step to result

Fast

- 10 min for drug measurement
- 12 min for ADA^b measurement

Flexible

- Whole blood (Capillary / Venous)**
- Serum
- Plasma

Reliable

- Excellent correlation with i-Tracker² and LISA TRACKER³ reference kits
- Calibrated on NIBSC/WHO international standards**
- Validated on princeps and biosimilars

* Therapeutic Drug Monitoring

** depending on test: refer to the instruction for use



Time Resolved Fluorescence technology


Connectable to the LIS

Friendly user interface


Compact and lightweight

Simple and quick test procedure

Serum / Plasma



Whole blood



MEASUREMENT RANGE

Drug

Infliximab 0,2 – 50 µg/mL

Adalimumab 0,2 – 50 µg/mL

Golimumab 0,2 – 30 µg/mL

Anti-Drug Antibodies

Anti-Infliximab 4 – 250 AU/mL

Anti-Adalimumab 3 – 200 AU/mL

Anti-Golimumab 2 – 40 ng/mL

Anti-Infliximab Total Ab 8 – 250 AU/mL

ez-Tracker

Reference

ETx 002-24

ETx 003-24

ETI 003T-24

ETx 002-C / ETx 003-C

ETx 002-CAL / ETx 003-CAL

Designation

ez-Tracker Drug

ez-Tracker Anti-Drug Antibodies

ez-Tracker Anti-Infliximab Total Ab

ez-Tracker Drug / ADAb Controls

ez-Tracker Drug / ADAb Calibrators

Packaging

24 tests

24 tests

24 tests

2 x 1 mL

2 x 1 mL

x = Infliximab / Adalimumab / Golimumab / Vedolizumab / Ustekinumab / Etanercept*

*Available soon

