

IMMUNO-TROL

AUTOIMMUNITY



A complete range of ready to use, internal autoimmune Quality Control sera, multiparametric, compatible with several technologies (multiplex FIDIS™, ELISA, DOT,...) and CE marked



- Celiac Disease
- Vasculitis
- Connective Tissue Diseases
- Rheumatoid Arthritis
- Antiphospholipid Syndrome
- Autoimmune Thyroiditis
- Inflammatory Bowel Disease
- Pernicious Anemia
- Diabetes Mellitus Type I
- Autoimmune Liver Diseases

IMMUNO-TROL sera are multiparametric Internal Quality Controls (IQC) for the routine validation of Autoimmune antibody assay techniques including FIDIS™ (xMAP® technology), ELISA and DOT.

ASSIGNED RANGES

- Liquid control, ready to use
- Standardised controls
- Identical use and reactivity as patients samples

CONVENIENT

- Universal multiparametric controls
- Stable material suitable for intra and inter assay validation

ADVANTAGES

- Monitor potential drift of analytical systems
- Quality assurance for assay validation

A LARGE QUALITY CONTROL RANGE OF 41 PARAMETERS IN AUTOIMMUNITY

	dsDNA	SS-A 60 Kda	TRIM 21 (SS-A 52 Kda)	SS-B	Sm	Sm/RNP	U1-RNP	Scl-70	Centromere protein B	Jo-1	Ribosomal proteins	Cardiolipin IgG	β2 GPI IgG	Prothrombin IgG	Phospholipids IgG	Cardiolipin IgM	β2 GPI IgM	Prothrombin IgM	Phospholipids IgM	Animal and human RF IgM	CCP IgG	TPO	TG	R-TSH	Mitochondria	MPO	PR3	GBM	PCA	IF	Transglutaminase IgG	Glialin IgG	Deamidated Gliadin IgG	Transglutaminase IgA	Glialin IgA	Deamidated Gliadin IgA	ASCA IgA	ASCA IgG	GAD	IA-2	ZnT8			
CONNECTIVE 1																																												
CONNECTIVE 2																																												
APS IgG																																												
APS IgM																																												
RHEUMA																																												
THYRO																																												
CYTO																																												
VASCULITIS																																												
GASTRITIS																																												
CELIAC G																																												
CELIAC A																																												
ASCA																																												
DIABETES																																												

ORDERING INFORMATIONS

REFERENCE	DESIGNATION	PACKAGING
HM 059-2X1	IMMUNO-TROL NEGATIVE	2X1 mL
HM 060-2X05	IMMUNO-TROL CONNECTIVE 1	2X0,5 mL
HM 061-2X05	IMMUNO-TROL CONNECTIVE 2	2X0,5 mL
HM 062-2X05	IMMUNO-TROL APS IgG	2X0,5 mL
HM 063-2X05	IMMUNO-TROL APS IgM	2X0,5 mL
HM 064-2X05	IMMUNO-TROL RHEUMA	2X0,5 mL
HM 065-2X05	IMMUNO-TROL THYRO	2X0,5 mL
HM 066-2X05	IMMUNO-TROL CYTO	2X0,5 mL
HM 067-2X05	IMMUNO-TROL VASCULITIS	2X0,5 mL
HM 068-2X05	IMMUNO-TROL GASTRITIS	2X0,5 mL
HM 069-2X05	IMMUNO-TROL CELIAC G	2X0,5 mL
HM 070-2X05	IMMUNO-TROL CELIAC A	2X0,5 mL
HM 071-2X05	IMMUNO-TROL ASCA	2X0,5 mL
HM 072-2X1	IMMUNO-TROL DIABETES	2X1 mL



Manufacturer: BIOSYNEX (France)

In vitro diagnostic medical devices.
Carefully read the instructions for use.
Pictures may differ from actual products.

BIOSYNEX

Theradiag

BIOSYNEX S.A.

14 rue Amboise Croizat

77183 Croissy-Beaubourg • France

Tel. : +33 (0) 1 64 62 10 12

Fax : +33 (0) 1 64 62 09 66

E-mail : service.clients@biosynex.com

www.theradiag.com

BIOSYNEX

BIOSYNEX S.A.

22 Boulevard Sébastien Brant

67400 Illkirch-Graffenstaden • France

Tel. : +33 (0) 3 88 78 80 88

Fax : +33 (0) 3 88 78 76 78

E-mail : service.clients@biosynex.com

www.biosynex.com